# MoVeS

# OXYCYCLE 1

## User Manual

English .....	2
Nederlands .....	8
Deutsch .....	14
Français .....	20
Espagnol .....	26
Svenska .....	32
Suomi .....	38

AROV1

The specifications of this product may vary from this photo, subject to change without notice.

## Table of Contents

---

Safety Precautions .....	2
Overview Drawing .....	3
Parts List .....	4
Assembly Instructions .....	5
Operating the Controller .....	6
Operating the Resistance Control Knob .....	6
Workout .....	7

## Safety Precautions

---

Read all instructions carefully before operating. Retain this owner's manual and keep all original transaction receipts for future reference.


- Consult your physician before beginning your exercise program.
- If you experience any irregular physical conditions such as dizziness, severe muscle or joint pain or pain in the chest, stop exercising and consult your physician immediately.
- Keep children and pets away from the machine while in use.
- Do not wear loose clothing while exercising.
- Keep a minimum safety clearance during operation.
- Do not operate the exercise equipment if it is damaged.
- Set up Pedal Exerciser as directed in Assembly Instructions and place unit on a flat, stable, non-slip surface.
- Perform 3 to 5 minutes of warm-up and stretching before beginning your exercise.
- Each exercise should be performed in a controlled manner. Always start exercising slowly.
- Never exercise to exhaustion.
- Do not stand on the Pedal Exerciser.
- The Pedal Exerciser is designed as a home use exercise device.
- Users must have full control of their muscles to operate this unit. It may not be suitable for people who are quadriplegic or paraplegic. Please consult your doctor or physician.
- Exercise equipment is to be regularly inspected for maintenance for loose or broken parts.
- Have the unit serviced by a qualified service technician. Do not attempt to service the unit yourself.

### Storage Precautions:


Do not store the unit in any place where it will be subject to high or low temperatures. Do not expose to direct sunlight for extended periods of time. Do not store the unit in a humid or dusty area.

# Overview Drawing

## Parts List

(18) Hexagon Socket  
Button Head Bolt  
6 PCS

(19) Spring Washer  
6 PCS

Wrench  
1 PC




Allen Wrench  
1 PC

Part #	Description	Q'ty	Part #	Description	Q'ty
001	Main frame	1	025	Eyebolt M6x43	2
002	Rear Stabilizer Ø50x1.5Tx380	1	026	Hexagon Nut M10x8.5	2
003	Front Stabilizer Ø50x1.5Tx220	1	027	Tension Bracket	2
004	Magnet Bracket	1	028	Hexagon Nut M6x6	2
005	Tension Control Knob Bracket	1	029	Hexagon Nut M8x7.7	1
006	Belt Pulley Ø120x152	1	030	Washer Ø8.5x16x1.5T	2
007	Magnet 15x10x5T	7	031	Hexagon Socket Button Head Bolt M8x42	1
008	Spacer Ø14xØ10.5x4	2	032	Cross Recessed Truss Head Tap- ping Screw ST4.2x15	4
009	Flywheel	1	033	Cross Recessed Pan Head Tapping Screw ST4.2x15	7
010	Stabilizer End Cap	4	034	Single Sided Adhesive Tape	4
011	Square Sensor (L=300mm)	1	035	Single Sided EVA Adhesive Tape	2
012	Big Washer Ø75x2.0T	2	036	Cross Recessed Pan Head Tapping Screw ST2.9x10	2
013	Bearing Bracket Ø74xØ41.5x12.5	2	037	Spring Ø9.4xØ1.0x50	1
014	Right Crank 5" (9/16")	1	038	Belt 220PJ3	1
015	Left Crank 5" (9/16")	1	039	Crank Cap	2
016	Digital Display Computer	1	040	Right Cover	1
017	Hexagon Flange Bolt M8x20	2	041	Left Cover	1
018	Hexagon Socket Button Head Bolt M8x12	6	042	Right Hand grip 9/16"	1
019	Spring Washer Ø8.4xØ13.5x2.5T	6	043	Left Hand grip 9/16"	1
020	C Ring Ø17xØ15.7x1.0T	4	044	Right Foot Pedal (JD-5B)	1
021	Bearing 6003Z	2	045	Left Foot Pedal (JD-5B)	1
022	Tension Control Knob	1	046	Round Magnet Ø12x5.0t	1
023	Flywheel Axle Ø10x75	1	047	Non-Slip Mat	1
024	Cross Recessed Truss Head Bolt M5x12	10			

## Assembly Instructions

---

1. Hold the unit in an upright position. The controller side facing up. Insert the front stabilizer (3) into the main frame (1).
2. Attach the front stabilizer (3) to the main frame (1) with three bolts (18) and three spring washers (19).
3. Insert the rear stabilizer (2) into the main frame (1).
4. Attach the rear stabilizer (2) to the main frame (1) with three bolts (18) and three spring washers (19).
5. Connect the left foot pedal (45) to the left crank (15) assembly. Thread it into the crank assembly in a counter-clockwise direction with the tool provided. (See figures 5 and 6). Connect the right foot pedal (44) to the right crank (14) assembly by threading it in a clockwise direction. Note: The left/right pedals and left/right cranks are marked with "L" and "R" (left and right)

## Operating the Digital Display Computer

---

### Computer Button Functions

Push the button to turn on the display.

Push the button to select a function.

Push and hold the button for 3 seconds to reset all functions except for the ODO function.

### Computer Display Functions & Measuring Range

SCAN		Automatically displays each function every 6 seconds.
TMR	[0:00~99:59]	Displays the exercise time.
KM	[0.0~199.9 km/h]	Displays the current training speed in kilometers / hour.
DIS	[0.0-999.9 km]	Displays the accumulative distance traveled during exercise.
CAL	[0.0-999.9]	Displays the calories burned during exercise.
ODO	[0.0-9999 km]	Displays the total accumulative distance traveled.
RPM	[0~9999 rpm]	Displays the current revolutions per minute.
STOP		The unit is stopped.
AUTO START		Start motion or push button.
DISPLAY OFF		After approximately 4 minutes without operation.

### Battery Replacement

Remove the display by using a slotted screwdriver.

Remove the battery with a slotted screwdriver. Replace with a new 3V lithium button cell battery.

## Operating the Resistance Control Knob

---

### 8-level resistance control knob

The resistance can be adjusted by turning the resistance control knob.

To increase the resistance, turn the resistance control knob in a clockwise direction. To decrease the resistance, turn the resistance control knob in a counterclockwise direction.

### How to use your OxyCycle 1 Pedal Exerciser.

---

When using the OxyCycle 1, the proper way to exercise is to set the desired resistance by turning the resistance control knob, and then pedaling to get the best results. Pedaling at excessive speeds or applying excessive force will damage the product or may cause injury to persons. Please start exercising at a slow pedaling speed.

## Workout


---

### Leg Exercise

Place the Pedal Exerciser on a non-slip mat and sit comfortably in a chair, the device directly in front of you. The knee joint should not be flexed less than 90 degrees. Wear rubber soled shoes and adjust foot pedal straps to the proper length. Barefooted use of the OxyCycle is not recommended.

### WARNING

Do not stand on the unit to exercise. Only exercise in a seated position. Make sure the seat does not move or pivot. Do not sit in a chair with casters.


NOTE: Use the foot pedals for leg exercises.

### Arm Exercise

Place the unit directly in front of you on a table top (on the included anti-slip mat). Use the hand grips for arm exercises.

[WWW.MOVES-YOU.COM](http://WWW.MOVES-YOU.COM)

MVS In Motion  
Westdijk 150  
2830 Tisselt (Willebroek)  
Belgium  
[www.mvs-in-motion.com](http://www.mvs-in-motion.com)

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© 2021, MVS In Motion. All rights reserved.