

Mobiderm Autofit lower limb

Main product features

January 2017

NEW IN 2017: Launch of the lower limb version

- > Range: 2 models proposed
 - Sock (knee high)
 - Thigh high stocking (full leg)
 - Standard equivalence to the existing made-to-measure Mobiderm models

> Technical choices

- Product adapted both to right/left limb => item codes optimisation
- 4 lengths (vs. 2 for the existing upper limb) => to match a larger patient size range and a mixt population (M/W concerned by lower limb LO)
- 4 sizes (acc. to circumferences), each one with a wider fitting possibility => to best answer the morphology variation of each patient
- New velcro quality & associated velvet material => reinforced hold!



NEW IN 2017: Launch of the lower limb version

Leg part: Mobiderm big

blocks

Foot part : Mobiderm small

blocks

Semi open garment: Easy to put on, easy to take off



Fitting marks: precision and control

Anti slip material: to avoid slipping on the floor



New Velcro straps & velvet material : reinforced hold!

Double raw of straps (to fold forward): wider fitting range/ avoids discomfort of the straps along the controlateral limb

Size range: Sock

16 item codes




4 sizes

	Circonférence de cheville Ankle level (cB)	Circonférence de mollet Calf level (cC)	Circonférence sous le genou Below knee level (cD)
T1	16 - 24 cm	24 - 35 cm	25 - 36 cm
T2	20 - 28 cm	32 - 43 cm	32 - 43 cm
T3	24 - 32 cm	40 - 51 cm	39 - 50 cm
T4	28 - 36 cm	48 - 59 cm	46 - 57 cm

Each size is proposed in 4 lengths

	Longueur mi-jambe Half leg length (/ D)
Court / Short	38 - 40 cm
Normal / Normal	40 - 42 cm
Long / Long	42 - 44 cm
Très long / Xtra long	44 - 46 cm



(CD) Floor-to-knee height (measurement taken below the bend of the knee, back side)

Size range: Thigh high stocking

16 item codes




4 sizes

	Circonférence de cheville Ankle level (cB)	Circonférence de mollet calf level (cC)	Circonférence de cuisse Thigh level (cG)
T1	16 - 24 cm	24 - 35 cm	46 - 57 cm
T2	20 - 28 cm	32 - 43 cm	53 - 65 cm
Т3	24 - 32 cm	40 - 51 cm	61 - 72 cm
T4	28 - 36 cm	48 - 59 cm	68 - 80 cm

Each size is proposed in 4 lengths

	Longueur de jambe Full leg length (/CG)
COURT / short	67 - 72 cm
Normal / Normal	72 - 77 cm
LONG / Long	77 - 82 cm
Très long / xtralong	82 - 87 cm



(&G) Floor-to-thigh height (at the gluteal fold)