



# micro 2<sup>m</sup>



## micro 72<sup>™</sup> Manual Muscle Tester

- The original digital dynamometer
- Ergonomic shape and lightweight design makes the microFet 2 easy to use
- · Compact design fits easily in the hand
- The industry standard in durability, reliability and accuracy
- Used by NASA, OSHA, and leading universities and research institutions worldwide
- Recognized by medico-legal community for documentation, accuracy and validity
- Versatile attachments ensure proper muscle isolation for increased testing validity
- Easy to read LCD display located directly on gauge shows peak force and duration
- For use as standalone or with our Clinical or FET Data Collection Software



EVALUATION TOOLS TO MEASURE, OBJECTIFY AND DOCUMENT HUMAN PERFORMANCE AND CARE



MT/ ROM



Pinch & Grip



Intelligent ROM

**Dual ROM** 





#### micro 2 TM





### micro 72 Specifications

- Patent # 5090421
- FDA Registered
- Size: 4" x 1 <sup>3</sup>⁄<sub>4</sub>"
- Weight: 0.8 lbs.
- Accurate to 1% full scale
- · Selectable units of measure: Pounds, Newtons or KGF.
- Two threshold settings for MT
  - Low threshold- 0.8 lbs. to 300 lbs. in 0.1 lb. increments
  - High threshold- 3.0 lbs. to 300 lbs. in 0.1 lb. increments
- Self-activating sleep mode to extend battery life
- 1 year warranty (Extended Warranties Available)

#### micro**FE72**<sup>™</sup> **Product Includes:**

- microFET2
- 3 attachments: Large, Digit and Flat
- Optional attachments available
- User Manual

- Upper and Lower Body Test Recording Pads

- Carrying Case
- Calibration Certificate

Muscle Test Wallchart



8020 South 1300 West West Jordan, UT 84088 USA

microFET products are manufactured to the highest standards of accuracy, reliability, and durability



VISIT US!

(800) 678-7888 or (801) 572-6500 contact@hogganhealth.com