


TD-1242

Non-Contact Forehead Thermometer



**Dual-Color
Backlit Display**



3 measurement modes

- ▶ Forehead General
- ▶ Forehead Child
- ▶ Surface mode



**Measurement
Distance 5 cm**

SPECIFICATION

Measuring Mode	Forehead (General/Child)	Surface
Measurement Range	32°C to 43°C (89.6°F to 109.4°F)	0°C ~ 100°C (32°F ~ 199°F)
Accuracy	35°C ~ 42°C (±0.2°C) ; < 35°C or > 42°C (± 0.3°C)	±1°C
Hi/Lo Indicator	32°C 38°C/37.7°C 43°C Lo Normal Fever Hi	0°C 100°C Lo Normal Hi
Operating Condition	10°C to 40°C (50°F to 104°F), 30-85% RH (non-condensing)	
Storage / Transportation Condition	-20°C to 60°C (-4°F to 140°F), 30-85% RH (non-condensing)	
Measuring Distance	5 ± 2 cm	
Resolution	0.1°C / 0.1°F	
Measurement Units	°C (preset) / °F	
Power Source	2 x 1.5V AAA alkaline battery	
Battery Life	> 3000 measurements	
Memory Capacity	30 measurements	
Power Saving	30 seconds	
Dimension	155.46mm x 40.14mm x 39.45mm	
Weight	61.8g without batteries	

SHIPPING INFO

Units per Carton	40
Gross Weight per Carton	9.2 kg
Carton size	39 (L) x 32.5 (W) x 36.5 (H) cm



TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288

www.taidoc.com