



CMA Testing and Certification Laboratories

廠商會檢定中心

TEST REPORT

Report No : AZ0014869(6)

Date : 09 Apr 2020

Application No : LZ007082(6)

Applicant :

Sample Description : One (1) submitted sample(s) stated to be :
Item Name : CHILDREN FACE MASK
Item No. : AN-002

Date Received : 30 Mar 2020.

Test Period : 30 Mar 2020 to 08 Apr 2020.

Test Requested	<u>Test Item</u>	<u>Result</u>
	1. Bacterial Filtration Efficiency (BFE)	Comply
	2. Particle Filtration Efficiency (PFE)	Comply

Test Result : Refer to the result pages for details

For and on behalf of
CMA Industrial Development Foundation Limited

Authorized Signature : _____ Page 1 of 3

Wan Leong Hang
Deputy Manager



CMA Testing and Certification Laboratories

廠商會檢定中心

TEST REPORT

Report No : AZ0014869(6)

Date : 09 Apr 2020

Application No : LZ007082(6)

Test Result :

No.	Test Item	Test Method	Limit	Test Result	Evaluation
1	Bacterial Filtration Efficiency (BFE)	YY 0469-2011, 5.6.1 (See Remark 2)	$\geq 95\%$	99.9 %	Comply
2	Particle Filtration Efficiency (PFE)	YY 0469-2011, 5.6.2 (See Remark 3)	$\geq 30\%$	95.3 %	Comply

Remark :

- (1) Notation :
% denotes Percentage
- (2) Bacterial filtration efficiency requirement stated in YY 0469-2011, 4.6.1 : $\geq 95\%$
- (3) Non-oily particulate filtration efficiency requirement stated in YY 0469-2011, 4.6.2 : $\geq 30\%$



CMA Testing and Certification Laboratories

廠商會檢定中心

TEST REPORT

Report No : AZ0014869(6)

Date : 09 Apr 2020

Application No : LZ007082(6)

Appendix



***** End of Report *****