



CMA Testing and Certification Laboratories

廠商會檢定中心

TEST REPORT

Report No. : AZ0015777(5)

Date: 20 Apr 2020

Application No. : LZ007892(4)

Applicant :

Sample Description : Submitted sample(s) stated to be :
Item Name : Premium Face Mask
Item No. : AN-001

Date Received : 09 Apr 2020.

Test Period : 09 Apr 2020 to 20 Apr 2020.

Test Requested : Test Item
1. Bacterial Filtration Efficiency (BFE)
2. Particle Filtration Efficiency (PFE)

Test Details :

No.	Test Item	Test Method	Limit	Test Result	Evaluation								
1	Bacterial Filtration Efficiency (BFE)	YY 0469-2011, 5.6.1 (See Remark 2)	$\geq 95\%$	<table border="0"> <tr> <td></td> <td style="text-align: center;">%</td> </tr> <tr> <td>1</td> <td style="text-align: center;">99.9</td> </tr> <tr> <td>2</td> <td style="text-align: center;">99.8</td> </tr> <tr> <td>3</td> <td style="text-align: center;">99.9</td> </tr> </table>		%	1	99.9	2	99.8	3	99.9	Comply
	%												
1	99.9												
2	99.8												
3	99.9												
2	Particle Filtration Efficiency (PFE)	YY 0469-2011, 5.6.2 (See Remark 3)	$\geq 30\%$	95.2%	Comply								

- Remark : (1) Notation :
 % denotes Percentage
(2) Bacterial filtration efficiency requirement stated in YY 0469-2011, 4.6.1 : $\geq 95\%$
(3) Non-oily particulate filtration efficiency requirement stated in YY 0469-2011, 4.6.2 : $\geq 30\%$

For and on behalf of
CMA Industrial Development Foundation Limited

Authorized Signature : _____

Tin Chun Shing
Assistant Manager

Page 1 of 2



CMA Testing and Certification Laboratories

廠商會檢定中心

TEST REPORT

Report No. : AZ0015777(5)

Date: 20 Apr 2020

Application No. : LZ007892(4)

Appendix



***** End of Report *****