



NEW RevoFit2™

With the new **RevoFit2™**, you can build in user adjustability *anywhere* on a socket.

Micro/Macro Adjustments

Each “click” of the dial provides 1mm of movement for targeted compression. When used correctly, the system can adjust for up to 10 ply of sock changes.

Adjustable Through Clothing

Operate the dial through clothing without removing or adding layers to adjust fit of the socket.



Rapid Donning & Doffing

Releasing the tension in the system makes donning and doffing easy for bulbous ends or large volume fluctuations.

Seated Comfort

Easily release the dial through clothing to make sitting more comfortable.

Need user-friendly suspension solutions? Order RevoLock™ Lanyard (PK3000-220-05).

The NEW RevoFit2™ design combines RevoFit™ Direct and RevoFit™ Versa into a single lamination dummy. RevoFit2™ kits are available for both diagnostic check sockets and definitive wet lamination.

RevoFit2™ Diagnostic Kit: SKU #PK2150-320-05



RevoFit2™ Lamination Kit: SKU #PK2000-320-05



Lamination 10-Pack: SKU #PK2010-320-05



Save money and have the right kit in stock by ordering a 10-pack of RevoFit2™ Lamination Kits.

For additional tech specs, please go to: www.ClickMedical.co/RevoFit2

Replacement Parts

For additional tech specs, please go to:
www.ClickMedical.co/RevoFit2

A full line of replacement parts allows the system to remain functioning like new throughout the lifetime of the socket.

SKU	Item No.	Details
RevoFit™ Replacement Lace Pack 	PK0060-000-10	Includes: <ul style="list-style-type: none"> • 7ft (2.1m) of 500-lb RevoFit™ Lace • Lace Feeder Tool • T6 Tool <p>We recommend replacing lace every 6 months or at any sign of abrasion to maintain proper system function.</p>
RevoFit™ Replacement Dial 	RP0100-320-05	Includes: <ul style="list-style-type: none"> • RevoFit2™ Dial • 7ft (2.1m) of 500 lb RevoFit™ Lace • Lace Feeder Tool • T6 Tool <p>We recommend replacing lace at time of dial replacement.</p>
RevoFit™ Bulk Lace	LOK014	Includes: <ul style="list-style-type: none"> • 10m of RevoFit™ Lace