



With RevoLock™ Lanyard, you can build in user-friendly suspension within a socket.

Automatic Alignment

RevoLock™ Lanyard automatically aligns the limb in the socket, eliminating pin lock frustrations.

Donn while Seated

RevoLock™ Lanyard gives patients the ability to fully draw the limb into the socket, even when seated. The powerful Boa® dial offers 2:1 mechanical leverage, giving the user the ability to draw in the liner even with limited strength.





Adjustable through Clothing

The dial can be operated through clothing, allowing the patient to adjust suspension or limb position as their volume changes.

Reduced Pistoning

Adjustable tension allows for limb to be properly suspended at all times, reducing pistoning.



Need adjustable socket solutions?
Order **RevoFit2™ Diagnostic Kit** (PK2150-320-05) and **Lamination Kit** (PK2000-320-05).

SKU	Item No.	Details
RevoLock™ Lanyard Kit 	PK3000-220-05	<ul style="list-style-type: none">• Diagnostic or Definitive: Serves as both a diagnostic and definitive socket kit.• Lower Extremity Locking Liner: Kit comes with threaded insert for lower extremity locking liners.• Standard Fabrication: RevoLock™ Lanyard is straightforward to assemble, easy to fabricate, and simple to incorporate within a traditional carbon fiber lamination.• Flexible Mounting: The dial can be mounted anywhere to accommodate the patient's functional needs and is easily adjustable through clothes and during activity.
Upper Ex Threaded Insert (6mm) 	RP3003-000-05	<ul style="list-style-type: none">• Upper Extremity Compatible: M6 Threaded Insert is available as an accessory for use with upper extremity or pediatric locking liners.

For additional tech specs, please go to: www.ClickMedical.co/RevoLockLanyard

Replacement Parts

For additional tech specs, please go to:
www.ClickMedical.co/RevoLockLanyard

SKU	Item No.	Details
RevoLock™ Replacement Lace 	RP3001-000-05	Includes: <ul style="list-style-type: none"> • 4ft (1.2m) of 300-lb RevoLock™ Lace • Lace Feeder Tool • T6 Tool <p>We recommend replacing lace every 6 months or at any sign of abrasion to maintain proper system function.</p>
RevoLock™ Replacement Dial 	RP3000-220-05	Includes: <ul style="list-style-type: none"> • RevoFit2™ Dial • 4ft (1.2m) of 300-lb RevoLock™ Lace • Lace Feeder Tool • T6 Tool <p>Dial snaps into laminated base plate. We recommend replacing lace at time of dial replacement.</p>
RevoFit™ Bulk Lace	LOK020	Includes: <ul style="list-style-type: none"> • 10m of RevoFit™ Lace