

PRODUCT CATALOG

FOOT ORTHOPEDICS

2025

MORE THAN 85%
OF STOCK ITEMS
READY FOR
DISPATCH



BAUERFEIND – BY PEOPLE FOR PEOPLE. FOR A BETTER QUALITY OF LIFE.

Bauerfeind AG is one of the world's leading manufacturers of medical aids. With our products, we enable people to be more active and thus make an important contribution toward maintaining and restoring health.

MEDICAL KNOW-HOW AND TECHNICAL EXPERTISE

Bauerfeind combines decades of expertise in technology and manufacturing with innovation: the tradition of a medium-sized, family-run business with a global mindset. Our specialists have worked very closely with physicians, scientists and developers right from the start. The result: patented, award-winning products and services.

EFFECTIVE VARIETY

We have taken on the mission of using modern orthotic blanks to make leading an active and pain-free life a reality, even well into old age. The products have been tried-and-tested on the market for many years and they are characterized by unique material properties.

Orthopedic foot orthoses relieve and stabilize the feet, knees, hips, and back. They have a biomechanical effect, helping the foot to function naturally. Pain can therefore be alleviated, mobility improved, and the overall quality of life increased.

UNCOMPROMISING QUALITY

Bauerfeind products don't simply meet standards – they fulfill the most stringent requirements when it comes to quality, effectiveness and wearing comfort. This is confirmed by regular inspection procedures, for which our products are permitted to bear the seals and marks after successful testing. Our orthotic blanks are developed and manufactured in Germany.

FROM THE DEVELOPMENT PHASE RIGHT THROUGH TO INDIVIDUAL ADJUSTMENT

Medical aids can only be 100 % effective if they fit properly. As different feet need different types of support, Bauerfeind offers numerous different types of foot orthoses, to cater for every requirement.



* except GloboTec Junior



CONTENTS

THE BAUERFEIND SERVICE YOU CAN RELY ON IT

> Page 6

ACADEMY COMPETENT. GOAL-ORIENTED. SUCCESSFUL.

> Page 8

MEASUREMENT TECHNOLOGY QUALITY THAT CAN BE MEASURED – AND MADE TO MEASURE

> Page 10

THE TECHNOLOGIES ADAPTED TO SUIT THE ANATOMY AND SYMPTOMS OF THE FEET

> Page 12

PRODUCT OVERVIEW AT A GLANCE

> Page 14

SPECIAL TREATMENT PURPOSES FOR STANDARD SYMPTOMS

> Page 16

THE RIGHT FOOT ORTHOSIS FOR EVERY FOOT

SYNTHETIC ORTHOTIC BLANKS



- 18 ErgoPad® weightflex
- 18 ErgoPad® weightflex 2*
- 22 GloboTec® Comfort
- 26 ErgoPad® redux heel
- 26 ErgoPad® redux heel 2*
- 28 ErgoPad® redux hallux
- 30 ErgoPad® work
- 30 Secosol® dynamic
- 34 GloboTec® Patricia®
- 37 GloboTec® PP-Flex (Polyflex)
- 38 GloboTec® Starflex®
- 40 GloboTec® Junior*
- 41 GloboTec® Isoflex
- 42 GloboTec®
Plastic / leather (Active)
- 43 GloboTec® Globoplast
- 45 TRIactive® Select
- 45 TRIactive® Shape

* medical devices

All products (except specifically labeled items) are stock items and are dispatched on the same day if ordered by 3 p.m.

SOFT FOAM ORTHOTIC BLANKS



- 48** TRlactive® Comfort
- 48** TRlactive® Soft
- 49** ErgoPad® Soft style, complete+
- 50** ErgoPad® system:s
- 52** GloboTec® Soft

MILLED ORTHOTIC BLANKS



- 62** BODYTRONIC® ID:CAM
- 64** ErgoPad® ID:
- 67** Milled sports foot orthoses
- 68** ErgoPad® ID:CoreX
- 70** ErgoPad® ID:Sandal

CORK ORTHOTIC BLANKS



- 74** GloboTec® Elastra-Cork
- 78** GloboTec® Alpha-Cork
- 79** GloboTec® Flexible Cork

VISCOELASTIC FOOT ORTHOSES



- 82** ViscoHeel®
- 82** ViscoHeel® K
- 83** ViscoSpot®
- 84** ViscoBalance®
- 85** ViscoPed®
- 85** ViscoPed® S

MATERIALS AND ADDITIONAL FEATURES



- 88** PedoPrint® footprint device
- 89** Imprint foam
- 90** Metatarsal pads
- 92** Cushions
- 93** Wedges
- 94** Heel cushions
- 94** Heel spur cushions
- 95** Small pads
- 96** Rubber foot orthoses
- 97** Top cover materials
- 100** Plate materials

OTHER BAUERFEIND RANGES

> Page 106

THE BAUERFEIND SERVICE YOU CAN RELY ON IT

ONLINE SHOP SERVICE


When ordering from the Bauerfeind online shop, you can benefit from reduced shipping costs and take advantage of various features such as the order status, availability display or patient files. To create a customer account, register at:

experts.bauerfeind.com/register



Order online now:

experts.bauerfeind.com




YOUR BENEFITS AT A GLANCE

- + 24 h online availability
- + Accessible anywhere using a web browser
- + Intuitive order process
- + Up-to-date order status
- + Low shipping costs
- + Smart search function
- + Digital order process
- + Direct connection to Bauerfeind's BODYTRONIC measurement technology
- + Patient file to store measurement and treatment data
- + Support via customer service, e-mail, contact form:
experts.bauerfeind.com/contact

CUSTOMER SERVICE FOR FOOT ORTHOSES

Benefit from our expertise and get advice.

We are available by phone from 8.00 a.m. to 6.00 p.m. from Monday to Friday. You can reach us 24 / 7 by e-mail and fax.



T +49 (0) 800 0 0105 30
F +49 (0) 800 0 0105 35
E info@bauerfeind.com

SERVICE GUARANTEES FOR MEDICAL RETAILERS IN GERMANY



SHIPPING ON THE SAME BUSINESS DAY

If you place your order by 3.00 p.m., we will ship in-stock items on the same day (large orders excluded). We are happy to accommodate your express delivery requirements, too.



FREE TELEPHONE AND FAX HOTLINES

We are happy to help you!



100% SERVICE AVAILABILITY

As competent contact partners, we are available by phone from 8.00 a.m. to 6.00 p.m. from Monday to Friday. You can reach us around the clock by fax, e-mail, or via electronic order systems (except in the event of system faults and maintenance work).



ELECTRONIC ORDER SYSTEMS

Place your orders easily, quickly and with reduced shipping costs: Using wheel-it or the various functions available in our online shop, such as the order status and availability display.



EXPRESS MANUFACTURE FOR QUICK CARE

Urgent, individually modeled milled foot orthoses – without intarsia – that are ordered before 10:00 a.m., will be completed and shipped on the same business day.



PRECISE DELIVERY DATE

We will specify a precise delivery date for your orders.



NO SHIPPING COSTS

We do not charge any shipping costs for order values over € 200 or for custom-made solutions. You can include all items from any of the product ranges in a single order.



FREE GOODS COLLECTION

If you fax us a completed return form, we will come and collect the goods to be returned (subject to inspection) – at our own expense.



365 DAY RETURNS

In accordance with our General Terms and Conditions, we will accept returns of stock items in faultless condition within 365 days following delivery.



RETURNS PROCESSED WITHIN FIVE DAYS

We will process your returns within five working days of receiving the goods.

BAUERFEIND ACADEMY

COMPETENT. GOAL-ORIENTED. SUCCESSFUL.

Would you like customers to clearly
favor your medical retailer and
reward your team's know-how?

We will be happy to support you!

We offer qualified training with strong practical
relevance all over Germany, with personal attendance
or online, for the following categories:

MEDICAL KNOW-HOW AND TECHNICAL EXPERTISE

"Basic foot orthosis training"

MARKETING / SALES

"The successful sale of medical products"

SPECIALIST SUBJECTS


"Workshop on digital orthotic care with
Bodytronic ID:CAM"

HUMAN RESOURCES

"Specialist medical supply retailer for
career changers and young professionals"
"Trainee camp"


Learn in your own time,
independent of location,
with our online courses.

bfzr.de/2-online



For the complete training
program, see

bfzr.de/2-akademie



Further information and registration
forms are available from your sales
representative or the Bauerfeind
Academy:

T +49 (0) 800 001 05 01

F +49 (0) 800 001 05 02

E akademie@bauerfeind.com

BENEFITS AT A GLANCE

+ TRAINING WITH STRONG PRACTICAL RELEVANCE

The training sessions offered by our academy will always include medical information to improve understanding and ensure increased technical expertise. Well-qualified trainers, such as orthotists, physical therapist and lymph therapists, will provide professional guidance through the different subjects, always considering the relevant target group.

+ IN-HOUSE TRAINING SESSIONS ACCORDING TO YOUR SPECIFICATIONS

Are you interested in personalized training that is tailor-made to your specific requirements? You can book your training session according to your specifications. You choose the subject, the location and the time. Get in touch!

+ CERTIFICATION

Our training sessions are certified in accordance with IQZ and AGOS.

+ BODYTRONIC EXPERT ADVISOR – OUR SPECIAL TIP

Become a "Bodytronic EXPERT ADVISOR" and convince your customers with digital measurement systems for your everyday customer service. Learn how to integrate digital measurement technology into consultation and sales, effortlessly and with confidence. You can measure discerning customers to provide products with a perfect fit and impress them with sophisticated technology.



QUALITY THAT CAN BE MEASURED – AND MADE TO MEASURE.

BAUERFEIND MEASUREMENT TECHNOLOGY



TOUCHPOINT THE ATTRACTIVE CONSULTATION STATION

- For straight-forward visualization and optimum treatment comfort
- 2D Rothballer scanner
- Quick measurement of the sole of the foot
- Documentation for traditional provision of care
- The recorded data can be used for the BODYTRONIC ID:CAM web application to model milled orthoses



BODYTRONIC® 140 THE DIGITAL BLUEPRINT

- 2D Rothballer scanner
- Quick measurement of the sole of the foot
- Suitable for mobile use as well as in-store
- Documentation for traditional provision of care
- The recorded data can be used for the BODYTRONIC ID:CAM web application to model milled orthoses





BODYTRONIC® 610

FULL-BODY SCANNER TO SUPPLY PRODUCTS FOR ALL EXTREMITIES

- Contact-free measurements to supply medical compression stockings, knee and ankle supports, knee orthoses as well as products for the back, arm and shoulder without the need for special measurement clothing
- Three-dimensional measurement of foam imprints; the recorded data can be used for the BODYTRONIC ID:CAM web application to model milled orthoses



BODYTRONIC® 150

WIRELESS MEASUREMENT OF DYNAMIC FOOT PRESSURE

- Insole measurement system by Medilogic
- For measuring dynamic pressure distribution across the sole of the foot
- For identifying the heel-to-toe movement and the body's center of gravity in the shoe
- Complete freedom of movement during the analysis, thanks to wireless use



BODYTRONIC® 410

PORTABLE 3D SCANNER

- Contact-free measurement to supply foot orthoses, with no need for special measurement clothing
- Three-dimensional measurement of foam imprints
- The recorded data can be used for the BODYTRONIC ID:CAM web application to model milled orthoses



THE TECHNOLOGIES

ADAPTED TO SUIT THE ANATOMY AND SYMPTOMS OF THE FEET


WIDESPREAD RELIEF

This technology targets anatomical weak points in the case of chronic heel and ankle pain or calcaneal spurs. The cut out provides relief for the heel area. The fan design relaxes the irritated area of the plantar aponeurosis. The core supports the foot arches, thus taking the excessive strain off the tendon.



RELIEF CORE TECHNOLOGY

- 1 Anatomically contoured with thin tapered edges:** Ensures the MTP joints are positioned outside the core, Avoiding any disruption of the heel-to-toe movement.
- 2 Fan-shaped cuts:** Relieves pressure on the fascia and tendon structure in the sole of the foot.
- 3 Arch support:** Raises flattened foot arches, thereby removing one of the causes of pain.
- 4 Cut out:** Relieves according to pathophysiological characteristics.
- 5 Heel cup:** Guides the foot with each step and promotes a biomechanically correct gait.



YOUR PERSONAL TOOLBOX

Using the Bodytronic ID:CAM system solution, milled orthoses can be modeled to offer a perfect, individual fit.



BODYTRONIC® ID:CAM

+ PROFITABILITY


- Free web application
- No maintenance or operating costs
- No need for your own materials warehouse

+ QUALITY

- Made-to-measure products in only 48 hours
- Precisely milled orthoses with surface properties that require little modification
- A reliable service for your customers


+ EFFICIENCY

- The perfect solution for every team
- Complete, workload-independent flexibility
- Intuitive operation and user interface



Try now at
idcam.bauerfeind.com






TRlactive®
TECHNOLOGY

A CLEVER CONNECTION

The TRlactive Technology uses a vertical serration system to combine materials with varying degrees of density, thus providing ideal support during every phase of walking.

- 1 **Guidance zone** 45 Shore
- 2 **Support zone** 35 Shore
- 3 **Cushioning zone** 25 Shore
- 4 **Functional zone**



TRlactive® Select rigidus


TRANSMISSION OF FORCE WITH THE X FACTOR

The "weightflex-X" integrated in the orthotic core supports the internal dynamics and torsional stability of the feet. The side sections guide occurring forces over two pivot points. The exchange of force supports the natural heel-to-toe movement of the feet and corrects the gait.


**WEIGHTFLEX®
TECHNOLOGY**

- Supports the inward rotation** of the forefoot
- Pivot points** of the "weightflex-X"
- Supports the longitudinal arch**

- 1 Hinge joints make the orthotic core bendable, so it can be adapted to suit a variety of heel heights.
- 2 The edging is anatomically contoured, leaving the metatarsophalangeal (MTP) joints free. The edges are thin.
- 3 The base of the fifth metatarsal bone is lower to provide relief. This takes into account the biomechanical requirements of this part of the foot.



The weightflex orthotic core is available in "soft," "medium," and "strong" versions to suit each patient's individual foot conditions, therapeutic goals, and body weight.






SUITABLE FOR VARIOUS HEEL HEIGHTS

Whether you are wearing ballet flats or heels – thanks to the integral hinges in the forefoot and heel areas, the orthotic core easily adapts to any heel height up to 6 cm.



THE RIGHT FOOT ORTHOSIS FOR EVERY FOOT

ORTHOTIC BLANKS	SYNTHETIC MATERIAL					
	ErgoPad® weightflex	GloboTec® Comfort	ErgoPad® redux heel	ErgoPad® redux hallux	TRlactive® Select Line	TRlactive® Shape
						
SIDE	18 – 21	22 – 25	26 – 27	28 – 29	45	45
SHOE TYPE						
Business	●	●				
Sports		●				
Kids		●				
Comfort	●	●	●	●	●	
Leisure			●	●	●	●
SYMPTOMS						
Arch decreased / splay foot	●	●	●	●	●	○
Skew foot			○	○	●	●
High-arched foot (mild / moderate)	○	○	○	○	●	
Chronic heel and ankle pain		●*	●		●*	
Problems in the forefoot area					●*	
Symptoms of the MTP joint		●*		●	●*	
Achilles tendon problems					○	
Knee problems	○	○	○	○	●*	○
MODE OF ACTION / BENEFITS						
Supportive foot orthoses	●	●	●	●	●	●
Cushioning foot orthoses					●	
Soft cushioning foot orthoses					●	
Shell foot orthoses		●				●

● Preferred use

● Main benefits

○ Further possible applications

* depending on the structural design of the foot orthosis or version

	SOFT FOAM					MILLED
	TRIactive® Comfort Line	TRIactive® Soft	ErgoPad® Soft style, complete+	ErgoPad® system:s	GloboTec® Soft	ErgoPad® ID: 
PAGES	48	48	49	50	52–59	64–66
SHOE TYPE						
Business						
Sports	●	●			●	●*
Kids					●	●
Comfort	●	●	●	●	●	●*
Leisure	●	●	●	●	●	●
SYMPTOMS						
Arch decreased / splay foot	●	●	●	●*	●	●*
Skew foot	●	●	○	●*	●	●*
High-arched foot (mild / moderate)	●	●	○	●*	○	●*
Diab. foot syndrome					○	●*
Rheumatoid foot def.					○	●*
Chronic heel and ankle pain	●*		○	●*	○	○*
Problems in the forefoot area	●*	●	○	●*		
Symptoms of the MTP joint	●*					
Achilles tendon problems	○*			○	○	○*
Knee problems	○	○	○	○	○	○
MODE OF ACTION / BENEFITS						
Supportive foot orthoses						●
Cushioning foot orthoses	●		●	●	●	
Soft cushioning foot orthoses		●	●	●	●	●
Shell foot orthoses					●	●

● Preferred use

● Main benefits

○ Further possible applications

* depending on the structural design of the foot orthosis or version

SPECIFIC TREATMENT FOR STANDARD SYMPTOMS*

Symptoms:		What can be done:	We recommend:
DIABETIC FOOT SYNDROME Diabetic foot syndrome is a secondary disease in patients suffering from diabetes mellitus. As the condition progresses, sensitivity reduces increasingly, and they feel less pain in their feet. That is why they notice pressure points too late or not at all. Wounds can then develop that become painful and infected. If these are not treated, they may result in permanently open wounds.		Foot orthoses with soft cushioning can perfectly distribute across the entire sole of the foot, pressure that occurs in specific areas during walking and standing. This can protect toes and the balls of the feet from pressure points. This reduces the risk of pressure ulcers and injuries.	<ul style="list-style-type: none"> • ErgoPad ID:Diabetes
CHRONIC HEEL AND ANKLE PAIN Chronic heel and ankle pain is very uncomfortable and occurs as soon as load is put on the heel. The tendon insertion on the heel bone is often irritated or inflamed because of excessive strain. Excessive tension on the heel bone's plantar tendon can result in a thorn-like protrusion on the bone – the heel spur.	 <p>Graphic representation</p>	Specifically adapted orthopedic foot orthoses can help with chronic heel and ankle pain. They surround the heel and provide support while the painful area features a cut out and particularly soft cushioning. This provides relief and alleviates pain immediately. Potential inflammation can then be reduced in the long term.	<ul style="list-style-type: none"> • ErgoPad redux heel • TRIactive Select heel • TRIactive Comfort heel • GloboTec Comfort FS • ViscoHeel, ViscoHeel K • ViscoSpot • ViscoPed, ViscoPed S • GloboTec Comfort style heel • ErgoPad work:h • ErgoPad work:h+
OSTEOARTHRITIS OF THE KNEE Osteoarthritis of the knee is a slowly progressing process that manifests with knee pain when subjected to strain, such as taking the stairs, walking on uneven terrain, or carrying heavy loads. Pain behind the kneecap after extended periods of sitting, increased sensitivity in damp or cold weather, and crunching sounds when moving the knee are also typical.	 <p>Graphic representation</p>	In order to alleviate pain caused by strain and to allow activity, the leg axis can be corrected using specific foot orthoses, thus relieving the affected joints. To this end, the foot orthosis features either a raised inner or outer edge. This distributes load on the joint, meaning pain during activity can be reduced significantly. In addition, the foot is guided during the heel-to-toe movement.	<ul style="list-style-type: none"> • TRIactive Select OA • GloboTec Soft • ErgoPad:ID
HALLUX RIGIDUS AND HALLUX VALGUS If the metatarsophalangeal joint is partly or fully impaired in its mobility, the natural heel-to-toe movement of the foot will be affected and cause severe pain. A joint condition (hallux limitus), a misalignment (hallux valgus), or complete stiffening of the joint (hallux rigidus) are often the cause.	 <p>Graphic representation</p>	A supporting orthopedic foot orthosis with what is known as a rigidus spring removes pressure from the painful joint when the foot is put down and reinforces it at the same time to restore natural movement. The metatarso-phalangeal joint is effectively relieved and immobilized.	<ul style="list-style-type: none"> • ErgoPad redux hallux • GloboTec Comfort style hallux • ErgoPad:ID • TRIactive Select rigidus • TRIactive Comfort rigidus
RHEUMATISM The feet of patients suffering from rheumatism are often swollen or even stiff. Contact or movement hurts and the skin over the joints is very sensitive.		Soft-cushioning foot orthoses provide feet affected by rheumatism with a comfortably soft contact surface. Feet sensitive to pressure and impact are noticeably relieved. Additionally, soft-cushioning foot orthoses support the foot arches, thus promoting the natural heel-to-toe movement of the foot.	<ul style="list-style-type: none"> • GloboTec Soft

* Illustration for a better overview only.
 For further information, please see: bauerfeind.com

SYNTHETIC MATERIAL

SYNTHETIC ORTHOTIC BLANKS


THE EXPERTISE
LIES IN THE CORE



ErgoPad® weightflex

ErgoPad® weightflex 2*

Comprehensive information relating to weightflex technology can be found on page 13.



	ErgoPad® weightflex	ErgoPad® weightflex 2*
USES	<ul style="list-style-type: none"> • Stable support for mild to pronounced flat foot and splay foot as well as mild high-arched foot • Available for strong, moderate and gentle support of the foot • Suitable for business and leisure shoes with a heel height of up to 6 cm 	
FEATURES	<ul style="list-style-type: none"> • Cushioning layer 	<ul style="list-style-type: none"> • Coated with high-quality polyurethane soft foam
	<ul style="list-style-type: none"> • Microfiber top cover 	<ul style="list-style-type: none"> • Fabric top cover
	<ul style="list-style-type: none"> • Stimulates the muscles through sensorimotor spots in the PUR surface 	
INFORMATION	<ul style="list-style-type: none"> • Supports the foot arches • Lowering at the base of the fifth metatarsal bone • Cores in three different levels of density (soft, medium, strong) • Supports the natural torsional movement of the foot • Can be adapted to suit various heel heights by means of hinges in the front of the core area and at the heel • Low construction height • A choice of narrow or wide cores and foot orthosis contours 	
MODIFICATIONS	<ul style="list-style-type: none"> • To be ground on the sides (at least 2 mm distance from core) 	<ul style="list-style-type: none"> • Cutting along the marked line creates a ¾ or a short-sole version
	<ul style="list-style-type: none"> • Easy to glue (standard workshop adhesives) 	
	<ul style="list-style-type: none"> • Orthotic cores are thermo-formable at 80 to max. 100 °C (warming in a water bath is recommended) 	<ul style="list-style-type: none"> • Complete thermoplastic molding possible (oven set to 90 °C, approx. 2 minutes)
	<ul style="list-style-type: none"> • When modifying the orthotic core, we recommend that you apply a top cover (e.g. made from microfiber or leather) 	
CLEANING / CARE	<ul style="list-style-type: none"> • Surfaces can be wiped clean using a mild cleanser and a sponge • Not machine-washable • Do not expose to direct heat (heater, direct sunlight) • We recommend that the functional properties are checked after six months 	

Applies for both foot orthoses versions

* Please refer to the relevant instructions for use for all information.

ErgoPad® weightflex



SYNTHETIC
ORTHOTIC BLANKS

Cores, wide

VERSION	SIZE (**)	ITEM NUMBER
soft	40–47	3 18031 1 00000**
medium	40–47	3 18032 1 00000**
strong	40–47	3 18033 1 00000**

Cores, narrow

VERSION	SIZE (**)	ITEM NUMBER
soft	35–41	3 18041 1 00000**
medium	35–41	3 18042 1 00000**
strong	35–41	3 18043 1 00000**

Base layers, wide



TOP COVER	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
–	Long	medium	40–48	3 18532 5 78110**
–	Long	strong	40–48	3 18533 5 78110**

Base layers, narrow

TOP COVER	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
–	Long	soft	35–42	3 18541 5 78110**
–	Long	medium	35–42	3 18542 5 78110**

Soft, narrow, long



TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Microfiber anthracite	–	Latex pad	Long	–	35–42	3 18541 5 74520**
Microfiber anthracite	PUR cushion 1.5 mm, ¾-length	Latex pad	Long	–	35–42	3 18541 5 77430**
Microfiber anthracite	PUR cushion 1.5 mm, ¾-length	Loose latex pad	Long	Open at the front	35–42	3 18541 5 78425**
Microfiber anthracite	Globofit blue 1.5 mm (12° Shore A)	Latex pad	Long	–	35–42 ¹	3 18541 5 78412**
Microfiber blue	Globofit blue 2 mm (12° Shore A)	Latex pad	Long	–	35–42 ¹	3 18541 5 78166**

¹ Not a stock item. Please note that a longer delivery time will apply.

Medium, narrow, long

TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Calfskin, perforated	Globofit beige 2 mm (12° Shore A)	Latex pad	Short	–	35–42	3 18542 5 47210**
Sleek leather, firm, perforated 1.0–1.25 mm	–	Loose latex pad	Short	Open at the front	35–42	3 18542 5 57440**
Microfiber anthracite	PUR cushion 1.5 mm, ¾-length	Latex pad	Long	–	35–42	3 18542 5 77430**

Medium, wide, long

TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Calfskin, perforated	Globofit beige 2 mm (12° Shore A)	Latex pad	Short	–	40–48	3 18532 5 47210**
Sleek leather, firm, perforated 1.0–1.25 mm	–	Loose latex pad	Short	Open at the front	40–48	3 18532 5 57440**
Microfiber anthracite	PUR cushion 1.5 mm, ¾-length	Latex pad	Long	–	40–48	3 18532 5 77430**
Microfiber anthracite	Globofit blue 1.5 mm (12° Shore A)	Latex pad	Long	–	40–48 ¹	3 18532 5 78412**
Microfiber blue	Globofit blue 2 mm (12° Shore A)	Latex pad	Long	–	40–48 ¹	3 18532 5 78166**

Strong, wide, long

TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, firm, perforated 1.0–1.25 mm	–	Loose latex pad	Short	Open at the front	40–48	3 18533 5 57440**
Microfiber anthracite	PUR cushion 1.5 mm, ¾-length	Latex pad	Long	–	40–48	3 18533 5 77430**

ErgoPad® weightflex 2

Narrow, long



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Fabric frost green	PUR	√	Soft, with PUR foam coating	35–42	3 76541 5 00000**
Fabric frost green	PUR	√	Medium, with PUR foam coating	35–42	3 76542 5 00000**
Fabric frost green	PUR	√	Strong, with PUR foam coating	35–42	3 76543 5 00000**

wide, long



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Fabric anthracite	PUR	√	Soft, with PUR foam coating	40–48	3 76531 5 00000**
Fabric anthracite	PUR	√	Medium, with PUR foam coating	40–48	3 76532 5 00000**
Fabric anthracite	PUR	√	Strong, with PUR foam coating	40–48	3 76533 5 00000**

GloboTec® Comfort

GloboTec® Comfort, style

The smart solution

ITEM NUMBER

3 12350 5 87210**



USES

- Moderate support for flattened or slightly raised foot arches with limited space available inside the shoe
- Ideal for customers / patients who are using foot orthoses for the first time
- An efficient, modern product that comes in a range of different versions and is easy to adapt
- Suitable for business and leisure shoes

FEATURES

- Semi-elastic orthotic blank
- Flexible and soft edges
- Low construction height, available in long-sole and short-sole versions
- High level of wearing comfort thanks to moderate foot arch supports
- Selected models feature a stabilizing heel half-cup

INFORMATION

MODIFICATIONS


- Partially or completely thermo-formable (hot air blower or oven, 110 °C, approx. 20–30 sec. for partial thermo-forming, approx. 2 min. for complete thermo-forming)
- Easy to grind (grinding belt, grit 60–100, polishing brush on the blade disk)
- Easy to glue (standard workshop adhesives)

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge (except leather top covers)
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

GloboTec® Comfort, style THE DESIGN IN DETAIL

- 1 MICROFIBER TOP COVER**
provides a comfortable contact surface for the foot
- 2 CUSHIONING LAYER**
long-sole soft bedding in the "individual solution" and "soft solution" models
- 3 METATARSAL PAD**
supports the natural arch of the foot behind the metatarsophalangeal joints (not available for "individual solution" model)
- 4 FOREFOOT CUSHION**
provides soft bedding for sensitive feet (not available for "individual solution" model)
- 5 BUTTERFLY CORE**
supports the arch of the foot
- 6 HEEL CUSHION**
ensures a soft tread
- 7 BOTTOM COATING**



Regular longitudinal arch, narrow

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Non-woven	–	–	Long	35–46	3 12540 5 00000**
Calfskin perforated	–	✓	Long	35–46	3 12640 5 00000**
Microfiber beige	–	✓	Long	35–46	3 12640 5 01400**
Microfiber anthracite	–	✓	Long	35–46	3 12640 5 01403**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	✓	Long	35–46	3 12740 5 01400**
onSteam anthracite	Globofit beige 2 mm (12° Shore A)	✓	Long	35–46	3 12740 5 01422**

Regular longitudinal arch, medium width

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Non-woven	–	–	Short	35–48	3 12530 1 00000**
Microfiber blue	–	✓	Short	35–48	3 12630 1 01401**
Non-woven	–	–	Long	35–48	3 12530 5 00000**
Calfskin perforated	–	✓	Long	35–48	3 12630 5 00000**
Microfiber beige	–	✓	Long	35–48	3 12630 5 01400**
Microfiber anthracite	–	✓	Long	35–48	3 12630 5 01403**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	✓	Long	35–48	3 12730 5 01400**
onSteam anthracite	Globofit beige 2 mm (12° Shore A)	✓	Long	35–48	3 12730 5 01422**

Regular longitudinal arch, wide

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Non-woven	–	–	Long	35–48	3 12550 5 00000**
Microfiber beige	–	✓	Long	35–48	3 12650 5 01400**
Microfiber anthracite	–	✓	Long	35–48	3 12650 5 01403**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	✓	Long	35–48	3 12750 5 01400**
onSteam anthracite	Globofit beige 2 mm (12° Shore A)	✓	Long	35–48	3 12750 5 01422**

GloboTec® Comfort



Special versions, long

DESCRIPTION	TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sports	Microfiber blue	–	✓	Wide, Davos lasts	35–48	3 12180 5 01401**
Sports	–	TECtop blue perforated 2 mm (30° Shore A)	✓	Narrow, Davos lasts	35–42	3 12190 5 22210**
Sports	–	TECtop blue perforated 2 mm (30° Shore A)	✓	Wide, Davos lasts	35–48	3 12180 5 22210**
–	Non-woven	–	–	Narrow, high longitudinal arch, moderate heel cup	35–42	3 12590 5 00000**
–	Non-woven	–	–	Wide, high longitudinal arch, moderate heel cup	35–48	3 12580 5 00000**
–	Microfiber anthracite	–	✓	Wide, moderate heel cup	35–48	3 12680 5 01403**
–	Microfiber blue	–	–	High heel cup	20–44 ²	3 12890 5 01401**
Chronic heel and ankle pain	onSteam Anthracite	Globofit beige 3 mm (12° Shore A)	–	Wide, heel cushion	35–48	3 54750 5 85121**



GloboTec® Comfort, style

DESCRIPTION	TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
The elegant solution	Microfiber orange	Functional cushion	Latex pad	Narrow	35–42	3 12380 5 86622**
The elegant solution, extra narrow	Microfiber navy	Globofit light blue 1.5 mm + Functional cushion	Latex pad	Extra narrow	35–42	3 12350 5 50000**
The smart solution	Microfiber navy	Functional cushion	Latex pad	Medium width	38–48	3 12350 5 87210**
The soft solution	Microfiber navy	Globofit light blue 2 mm (12° Shore A) + Functional cushion	Latex pad	Medium width	35–48	3 12350 5 88010**
The individual solution	Microfiber navy	Globofit light blue 2 mm (12° Shore A)	Latex pad	Medium width, open at the front	35–48	3 12350 5 86015**
Hallux both sides	Microfiber navy	Globofit light blue 2 mm (12° Shore A)	Latex pad	Medium width	35–48	3 36883 5 00000**
Hallux left	Microfiber navy	Globofit light blue 2 mm (12° Shore A)	Latex pad	Medium width	35–48	3 36881 5 80000**
Hallux right	Microfiber navy	Globofit light blue 2 mm (12° Shore A)	Latex pad	Medium width	35–48	3 36882 5 80000**
Heel	Microfiber gray	TECfom light azure 3mm	Latex pad	Medium width	35–48	3 12350 5 86141**

² Double sizing

GloboTec® comfort sports*

DESCRIPTION	TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
walk	Microfiber navy	PowerWave walk 3 mm	Latex pad	–	35–48	3 13100 5 00000**
walk	Microfiber navy	PowerWave walk 3 mm	–	–	35–48	3 13100 5 00350**
play	Microfiber navy	PowerWave play 3 mm	Latex pad	–	35–48	3 13150 5 00000**
play	Microfiber navy	PowerWave play 3 mm	–	–	35–48	3 13150 5 00350**
run	Microfiber navy	PowerWave run 3 mm	Latex pad	–	35–48	3 13250 5 00000**
run	Microfiber navy	PowerWave run 3 mm	–	–	35–48	3 13250 5 00350**
run	Microfiber navy	PowerWave run 3 mm	–	Loose top layer	35–48	3 13250 5 22110**
golf	Microfiber navy	PowerWave golf 3 mm	Latex pad	–	35–48	3 13400 5 00000**

**Children's foot orthoses, long, high heel cup**

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Microfiber blue (animal prints)	–	–	Narrow heel	20–44 ²	3 12870 5 00000**
–	TECmoll blue 2 mm (35° Shore A)	–	Wide heel	20–44 ²	3 12880 5 82333**
Microfiber blue (animal prints)	TECmoll blue 2 mm (35° Shore A)	–	Wide heel	20–44 ²	3 12880 5 00000**
Microfiber anthracite	TECmoll blue 2 mm (35° Shore A)	–	Wide heel	20–44 ²	3 12880 5 01403**
Microfiber blue	–	–	–	20–34 ²	3 12890 5 01401**

² Double sizing

* Please refer to the relevant instructions for use for all information.

ErgoPad® redux heel

ErgoPad® redux heel 2*

Comprehensive information relating to Relief Core Technology can be found on page 12.




ErgoPad® redux heel

ITEM NUMBER
3 12450 5 00000**




	ErgoPad® redux heel	ErgoPad® redux heel 2*
USES	<ul style="list-style-type: none">• Suitable for standard and comfort shoe shapes• For the effective treatment of chronic heel and ankle pain, plantar fasciitis and plantar calcaneal spurs• For preventing chronic symptoms• For supporting flattened or slightly raised foot arches	
FEATURES	<ul style="list-style-type: none">• With high-quality multi-layer cushions• Available with a glued or unglued cushion layer	<ul style="list-style-type: none">• Coated with high-quality polyurethane soft foam• Requires little modification
	<ul style="list-style-type: none">• Low construction height• Cushioning in the painful, pressure-sensitive heel area• Synthetic orthotic core with hole notches and an intelligent fan design for relieving the plantar aponeurosis	
INFORMATION		
MODIFICATIONS	<ul style="list-style-type: none">• Easy to grind (grinding belt, grit 60–100, polishing brush on the blade disk)• Easy to glue (standard workshop adhesives)	<ul style="list-style-type: none">• Cutting along the marked line creates a ¾ or short-sole version
	<ul style="list-style-type: none">• The separate top layer allows the metatarsal pads to be freely positioned and makes modifications possible	
CLEANING / CARE	<ul style="list-style-type: none">• Thermo-formable (oven set to 110 °C, approx. 5 minutes)	
	<ul style="list-style-type: none">• Surfaces can be wiped clean using a mild cleanser and a sponge• Not machine-washable• Do not expose to direct heat (heater, direct sunlight)• We recommend that the functional properties are checked after six months	

 Applies for both foot orthoses versions

* Please refer to the relevant instructions for use for all information.

ErgoPad® redux heel THE DESIGN IN DETAIL

- 1 **MICROFIBER TOP COVER**
- 2 **CUSHION LAYER**
with PUR cushion
- 3 **METATARSAL PAD**
- 4 **SYNTHETIC CORE**
with cut-out and fan-shaped
plantar section
- 5 **BOTTOM COATING**
with PUR cushion



ErgoPad® redux heel

Glued versions



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Microfiber macchiato	TECmoll silver-gray 3 mm (35° Shore A)	Latex pad	Narrow	35–42	3 12440 5 00000**
Microfiber macchiato	TECmoll silver-gray 3 mm (35° Shore A)	Latex pad	Narrow, flat longitudinal arch	35–42	3 12441 5 00000**
Microfiber anthracite	TECmoll silver-gray 3 mm (35° Shore A)	Latex pad	Wide	38–48	3 12450 5 00000**
Microfiber anthracite	TECmoll silver-gray 3 mm (35° Shore A)	Latex pad	Wide, flat longitudinal arch	38–48	3 12451 5 00000**

Loose top layer



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Microfiber macchiato	TECmoll silver-gray 3 mm (35° Shore A)	Loose latex pad	Narrow, comfort shoe shape	35–42	3 12402 5 00440**
Microfiber anthracite	TECmoll silver-gray 3 mm (35° Shore A)	Loose latex pad	Wide, comfort shoe shape	38–48	3 12400 5 00440**

ErgoPad® redux heel 2



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Fabric frost green	PUR	✓	Narrow, with PUR foam coating	35–42	3 77540 5 00000**
Fabric anthracite	PUR	✓	Wide, with PUR foam coating	38–48	3 77530 5 00000**

ErgoPad® redux hallux

ITEM NUMBER

3 36942 5 00000**



USES

- Primarily suited for non-surgical (pre-operative and post-operative) treatment of hallux limitus and hallux rigidus (stiffening of the big toe)

FEATURES

- Designed to limit painful movements in the first metatarsophalangeal joint (MTP 1)
- A viscoelastic reinforcement element can be applied on one side (right or left, with the opposite side flexible) or on both sides
- The reinforcement element supports the entire forefoot (preventing forefoot misalignments)
- The ball of the big toe lies freely and is more deeply embedded
- Shock-absorbing PUR cushion layer (2 mm)
- Low construction height

INFORMATION

MODIFICATIONS


- Requires little modification
- The longitudinal arch support can be thermo-formed (hot air blower, 110 °C, 20–30 sec.)
- The pads can be freely positioned in versions that are open at the front

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

THE DESIGN IN DETAIL

- 1 TOP COVER AND PUR CUSHION LAYER**
High-quality microfiber top cover; shock-absorbing PUR cushion layer
- 2 METATARSAL PAD**
- 3 BUTTERFLY CORE**
optimal support for the metatarsal area without allowing the foot orthosis to become torsion-resistant
- 4 ORTHOTIC CORE**
- 5 REINFORCEMENT ELEMENT**
provides relief for the painful area; prevents supination misalignment of the forefoot; avoids levering out of the shoes
- 6 FRAME**
- 7 BOTTOM COATING**
Tried-and-tested grip film – non-slip, flexible layer in contact with the shoe



Glued versions, long


TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Microfiber beige	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on the left (L), narrow	35–42	3 36941 5 00000**
Microfiber brown	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on the left (L), wide	38–48	3 36931 5 00000**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on the right (R), narrow	35–42	3 36942 5 00000**
Microfiber brown	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on the right (R), wide	38–48	3 36932 5 00000**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on both sides (L / R), narrow	35–42	3 36943 5 00000**
Microfiber brown	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on both sides (L / R), wide	38–48	3 36933 5 00000**

**Open at the front, long**

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Microfiber beige	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on the left (L), narrow	35–42 ¹	3 36941 5 22992**
Microfiber brown	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on the left (L), wide	38–48 ¹	3 36931 5 22992**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on the right (R), narrow	35–42 ¹	3 36942 5 22992**
Microfiber brown	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on the right (R), wide	38–48 ¹	3 36932 5 22992**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on both sides (L / R), narrow	35–42	3 36943 5 22992**
Microfiber brown	Globofit beige 2 mm (12° Shore A)	Latex pad	Reinforced on both sides (L / R), wide	38–48	3 36933 5 22992**

¹ Not a stock item. Please note that a longer delivery time will apply.

ErgoPad® work orthotic blanks*



	ErgoPad® work:h ⁺	ErgoPad® work:h
USES	<ul style="list-style-type: none"> • Moderate support for flattened foot arches, in cases of foot, knee or back pain, if required • For (work-related) strain on the locomotor system while standing, walking, or carrying heavy loads, and with excessive physical strain or premature muscle fatigue • Suitable for work and safety shoes as well as leisure and comfort shoes 	
FEATURES	<ul style="list-style-type: none"> • Synthetic orthotic core with hole notches and an intelligent fan design for relieving the plantar aponeurosis • Adequate orthopedic support for the foot arches and the musculoskeletal system • Cushioning areas with durable elasticity under the heel and the lateral midfoot and forefoot area • Stabilizing heel cup for improved hindfoot stability • ESD-conductive orthotic blanks 	
INFORMATION		
MODIFICATIONS	<ul style="list-style-type: none"> • In accordance with the enclosed instructions for use • Easy to grind (grinding belt, grit 24–40, polishing brush on the blade disk) • Easy to glue (using a design-tested adhesive and construction materials) • Partially thermo-formable (hot air blower, 120 °C, max. 20 sec.) 	
CLEANING / CARE	<ul style="list-style-type: none"> • Surfaces can be wiped clean using a mild cleanser and a sponge • Not machine-washable • Do not expose to direct heat (heater, direct sunlight) • We recommend that the functional properties are checked after six months 	

Applies for both foot orthoses versions

LEGAL NOTE

The supply and use of foot orthoses for work safety footwear are subject to safety specifications in accordance with DIN EN ISO 20345 as well as DIN ISO 61340 because of possible ESD conductivity. Please check with the manufacturer of the work safety footwear whether the required type testing has been carried out. Bauerfeind cannot answer this question because shoe collections change and there are many manufacturers.

When manufacturing the foot orthoses for work safety footwear, the relevant instructions for use for processing the specific orthotic blanks must be complied with to prevent risks that may arise when the shoes are being worn. This compliance forms the basis for type testing performed by shoe manufacturers. Non-compliance with relevant specifications contained in the instructions for use results in the exclusion of any warranty claims associated with this.

* Please refer to the relevant instructions for use for all information.

THE DESIGN IN DETAIL

ErgoPad® work:h⁺

- 1 TOP COVER AND CUSHION LAYER**
Conductive PUR cushion layer with polyester top cover, ESD-conductive foam
- 2 TEAR-SHAPED METATARSAL PAD**
Supports the plantar arch behind the metatarsophalangeal joints

- 3 SYNTHETIC CORE – RELIEF CORE TECHNOLOGY**
Cut out in the heel area with heel pad for additional cushioning and plantar fan design to relieve the irritated plantar aponeurosis

- 4 BOTTOM COATING**
Tear-resistant and conductive bottom coating



THE DESIGN IN DETAIL

ErgoPad® work:h

- 1 TOP COVER AND CUSHION LAYER**
ESD conductivity based on carbon fibers woven into the polyester top cover

Soft cushion layer made of ESD-conductive foam (5 mm forefoot and heel cushioning)
- 2 ANATOMICALLY CONTOURED METATARSAL PAD**
supports the natural arch of the foot behind the metatarsophalangeal joints

- 3 SELF-SUPPORTING SYNTHETIC CORE**
Cut out in the heel area for additional cushioning; fan design on the sole to relieve the irritated tendon path (study results substantiate the high level of effectiveness of the orthotic core in terms of pain reduction – orthotic treatment for plantar fasciitis, study conducted by Prof. Dr. med. Markus Walther)*

- 4 BOTTOM COATING**
maximum, non-slip contact surface between the foot orthosis and the shoe; made of conductive EVA



* Walther, Markus et al.: Effect of different orthotic concepts as first line treatment of plantar fasciitis. In: Foot And Ankle Surgery, Vol. 19, Nr. 2, June 2013, S. 103-107

Secosol® dynamic

The certified foot orthosis
for safety footwear



Certified shoe
models:



The Secosol dynamic orthotic blanks have been optimized for use in work safety footwear. They relieve the feet during all standing activities. With a slight heel cup, a longitudinal arch support in the low and mid versions, and a metatarsal pad to relieve the forefoot, they support the feet and ensure secure footing. The orthotic blanks are thermo-formable and can therefore be individually adapted to the feet.

The certified shoe models are listed on the **Secosol shoe finder**.

Bauerfeind and
Hartmann Orthopädie+Sport GmbH
are retail partners.

SECOSOL® – EINE MARKE VON

Hartmann
Orthopädie+Sport

**ErgoPad® work:h+**

TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Fabric black	ESD foam	Latex pad	ESD foam	Long	35–48	3 11350 5 00000**
Fabric black	ESD foam	–	ESD foam	Long	35–48	3 11350 5 0035 0**
Fabric black	ESD foam	Latex pad	ESD foam	Flat longitudinal arch	35–48	3 11350 5 8000 0**
Fabric black	ESD foam	Latex pad	ESD foam	Loose top layer	35–48 ²	3 11350 5 0044 0**
Fabric black	ESD foam	Latex pad	ESD foam	Partially glued, long	35–48 ²	3 11350 5 0010 0**

**ErgoPad® work:h**

TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Fabric gray	ESD foam	Latex pad	ESD foam	Long	35–48	3 10350 5 00000**
Fabric gray	ESD foam	Latex pad	ESD foam	Flat longitudinal arch	35–48	3 10350 5 80000**
Fabric gray	ESD foam	–	ESD foam	Long	35–48	3 10350 5 00350**
Fabric gray	ESD foam	Loose latex pad	ESD foam	Long, open at the front	35–48 ²	3 10350 5 00100**
Fabric gray	ESD foam	Loose latex pad	ESD foam	Long, loose top layer	35–48 ²	3 10350 5 00440**

**ErgoPad® work:x**

TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Fabric gray	ESD foam	Latex pad	ESD foam	Medium, narrow	35–39 ²	3 10350 5 05820**
Fabric gray	ESD foam	Latex pad	ESD foam	Medium, narrow, loose top layer	35–39 ¹	3 10350 5 27128**
Fabric gray	ESD foam	Latex pad	ESD foam	Medium, wide	40–48 ²	3 10350 5 05820**
Fabric gray	ESD foam	Latex pad	ESD foam	Medium, wide, loose top layer	40–48 ¹	3 10350 5 27128**
Fabric gray	ESD foam	Latex pad	ESD foam	Strong, wide	40–48 ²	3 10350 5 05830**
Fabric gray	ESD foam	Latex pad	ESD foam	Strong, wide, loose top layer	40–48 ¹	3 10350 5 27138**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Available only in connection with the "Assumption of Liability Agreement".
Please contact our Foot Orthopedics Customer Service Team.



ErgoPad® work for safety shoes

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Velour black	ESD foam	√	Long	35–52	3 10600 5 00000**



Secosol® dynamic

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Fabric antistatic	EVA perforated	√	low	35–50	3 10500 5 80014**
Fabric antistatic	EVA perforated	√	mid	35–50	3 10500 5 00014**

 Manufacturer: Hartmann Orthopädie + Sport

GloboTec® Patricia®



GloboTec® Patricia®

Half-cup leather blank, long-sole

ITEM NUMBER

3 11062 5 00000**

USES

- For supporting flattened or slightly raised foot arches
- Suitable for people with a normal body weight
- Suitable for business and leisure shoes with limited space inside

FEATURES

- Flexible injection-molded orthotic core
- Thin materials, available in long-sole and short-sole versions
- Very thin at the heel and the edges
- Selected models feature a stabilizing heel half-cup

INFORMATION

MODIFICATIONS

- Partially thermo-formable Patricia® core (hot air blower, 130–150 °C, approx. 30 sec.)
- Easy to glue (standard workshop adhesives)

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge (except leather top covers)
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

Heel cup



TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 1.5mm	–	Latex pad	Clarino Nash	Narrow, short	35–44	3 11062 1 00000**
Sleek leather, soft 1.5mm	–	Latex pad	Clarino Nash	Narrow, long	35–42	3 11062 5 00000**
Microfiber blue	TECtop blue 3 mm (30° Shore A)	Latex pad	Carbon film	Narrow, long	35–42 ¹	3 11062 584239**
Sleek leather, soft 1.5mm	–	Latex pad	Clarino Nash	Wide, short	38–48	3 11063 1 00000**
Sleek leather, soft 1.5mm	–	Latex pad	Clarino Nash	Wide, long, square shape	38–48	3 11063 5 00000**
–	Lunasoft (four colors) 3 mm (40° Shore A)	Metatarsal pad from size 36	Non-woven	Narrow, long, loose supination wedges, heel cup	20–40 ²	3 11171 5 00000**
Microfiber blue	Lunasoft (four colors) 3 mm (40° Shore A)	Metatarsal pad from size 36	Non-woven	Narrow, long, loose supination wedges, heel cup	20–40 ¹⁺²	3 11171 5 85144**
Goatskin 0.8mm	TECmoll silver-gray 3 mm (35° Shore A)	Latex pad	Non-woven	Narrow, long	35–48	3 11562 5 00000**
Goatskin 0.8mm	TECmoll silver-gray 3 mm (35° Shore A)	Latex pad	Non-woven	Wide, long, square shape	38–48	3 11563 5 00000**

Sleek leather, soft 1.5 mm

METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Latex pad	Clarino Nash	Narrow, short (MTP joints free)	35–40	3 11052 1 00000**
Latex pad	Clarino Nash	Narrow, short (MTP joints free)	41–48	3 11053 1 00000**
Latex pad	Clarino Nash	Narrow, long	35–40	3 11212 5 00000**
Latex pad	Clarino Nash	Narrow, long	41–48	3 11213 5 00000**
Latex pad	Non-woven	Wide, long	35–40	3 11112 5 00000**
Latex pad	Non-woven	Wide, long	41–48	3 11113 5 00000**

Goatskin 0.8 mm + Globofit beige 3 mm (12° Shore A)

METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Latex pad	Non-woven	Narrow, long	34–40	3 11542 5 00000**
Latex pad	Non-woven	Narrow, long	41–48	3 11543 5 00000**
Latex pad	Non-woven	Wide, long	35–40	3 11282 5 00000**
Latex pad	Non-woven	Wide, long	41–48	3 11283 5 00000**

Microfiber beige + Globofit beige 3 mm (12° Shore A)

METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Latex pad	Non-woven	Narrow, long	35–40 ¹	3 11542 5 83297**
Latex pad	Non-woven	Narrow, long	41–48 ¹	3 11543 5 83297**

Goatskin 0.8 mm + TECsoft 3 mm (25° Shore A)

METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Latex pad	Non-woven	Narrow, long	35–40	3 11322 5 00000**
Latex pad	Non-woven	Wide, long	35–40	3 11332 5 00000**
Latex pad	Non-woven	Wide, long	41–48	3 11333 5 00000**

¹ Not a stock item. Please note that a longer delivery time will apply.² Double sizing

GloboTec® Patricia®



Goatskin 0.8 mm + TECmoll gray 3 mm (35° Shore A)

METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Latex pad	Non-woven	Narrow, long	34–40	3 11512 5 00000**
Latex pad	Non-woven	Narrow, long	41–46	3 11513 5 00000**
Latex pad	Non-woven	Wide, long	35–40	3 11532 5 00000**
Latex pad	Non-woven	Wide, long	41–48	3 11533 5 00000**



Microfiber anthracite + TECmoll gray 3 mm (35° Shore A)

METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Latex pad	Carbon film	Narrow, long	34–40 ¹	3 11512 5 84473**
Latex pad	Carbon film	Narrow, long	41–46 ¹	3 11513 5 84473**
–	Non-woven	Wide, long, open at the front	35–40 ¹	3 11532 5 35871**
–	Non-woven	Wide, long, open at the front	41–48 ¹	3 11533 5 35871**



Children's foot orthoses

TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather soft 1.5 mm	–	–	Clarino Nash	Narrow, long, heel cup	20–32 ² 33–34	3 11061 5 00000**
–	Lunasoft (four colors) 3 mm (40° Shore A)	–	Non-woven	Narrow, long, loose supination wedges, heel cup	20–34 ²	3 11171 5 00000**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

GloboTec® PP-Flex (Polyflex)



ITEM NUMBER

3 09339 1 03999**

USES

- For supporting flattened or slightly raised foot arches
- Suitable for people with a normal body weight

FEATURES

- Synthetic orthotic cores with various support functions
- Thin materials, available in long-sole and short-sole versions
- Thin at the heel and edges, some available with a stabilizing heel cup

INFORMATION

MODIFICATIONS

- The core is partially or completely thermo-formable (hot air blower or oven, 130 to 150 °C, 30 to 60 sec. for partial thermo-forming, approx. 3 to 4 min. for complete thermo-forming)
- Easy to grind (grinding belt, grit 60–100, polishing brush on the blade disk)
- Can be glued to other materials (standard workshop adhesives)

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	–	–	Wide, short, without top layer	24–34 ¹⁺² 35–48	3 09339 1 03999**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	Latex pad	Wide, long, open at the front	38–48	3 09931 5 01400**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

GloboTec® Starflex®



ITEM NUMBER:
3 05431 5 00000**

USES

- Stable support for flattened foot arches
- Suitable for average as well as higher body weights

FEATURES

- Stable injection-molded orthotic core
- Break-proof even under high loads
- Thin materials, available in long-sole and short-sole versions

INFORMATION

MODIFICATIONS

- Starflex Core is partially or completely thermo-formable (hot air blower or oven, max. 100 °C, approx. 60 sec. for partial thermo-forming, approx. 4 min. for complete thermo-forming)
- Easy to grind (grinding belt, grit 40–100, polishing brush on the blade disk)

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge (except leather top covers)
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

Cores

METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	Extra narrow, brushed	36–44 ²	3 05180 1 00000**
–	Narrow, brushed	34–46	3 05140 1 00000**
✓	Narrow	32–44	3 05141 1 00555**
–	Wide	24–34 ² 35–48	3 05130 1 00000**
✓	Extra wide, brushed	40–47	3 05111 1 00000**

Loose top layer

TOP COVER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather soft 1.5 mm	Loose latex pad	Short	Narrow, short	34–46	3 05341 1 00000**
Sleek leather soft 1.5 mm	Loose latex pad	Short	Narrow, long	34–46	3 05341 5 00000**
Sleek leather soft 1.5 mm	Loose latex pad	Short	Wide, short	24–34 ¹⁺² 35–48	3 05331 1 00000**
Sleek leather soft 1.5 mm	Loose latex pad	Short	Wide, long	24–34 ¹⁺² 35–48	3 05331 5 00000**
–	Latex pad	Short	Wide, short, without top layer	24–34 ¹⁺² 35–48	3 05339 1 00000**

Glued / partially glued versions

TOP COVER	CUSHION LAYER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft, perforated 1.25 mm	–	Latex pad	–	Narrow, long, glued	34–46	3 05441 5 00000**
Sleek leather, soft, perforated 1.25 mm	–	Latex pad	–	Wide, long, glued	24–34 ¹⁺² 35–48	3 05431 5 00000**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	Latex pad	Short	Narrow, long, open at the front	34–46	3 05941 5 01400**
Microfiber beige	Globofit beige 2 mm (12° Shore A)	Latex pad	Short	Wide, long, open at the front	24–34 ² 35–48	3 05931 5 01400**

¹ Not a stock item. Please note that a longer delivery time will apply.² Double sizing

GloboTec® Junior*



ITEM NUMBER:
3 10805 1 00000**

USES

- Correction and / or arch support for planovalgus foot
- Efficient, modern product that is easy to adapt
- Suitable for a wide range of children’s shoes

FEATURES

- Injection-molded orthotic core with cushion layer
- Thin materials, ¾-length version
- Very thin and flexible edges
- High, stabilizing heel cup
- Thin at the heel and edges

INFORMATION

MODIFICATIONS

- Requires little modification; usually only needs to be shaped to fit the shoes
- Partially thermo-formable (hot air blower, 120 °C, 20–30 sec.)

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Microfiber blue	TPE cushion	–	¾-length, high heel cup	22–40 ²	3 10805 1 00000**

* Please refer to the relevant instructions for use for all information.

² Double sizing

GloboTec® Isoflex



ITEM NUMBER
3 19631 5 00000**

USES

- For supporting flattened or slightly raised foot arches
- Suitable for people with a normal body weight

FEATURES

- Injection-molded orthotic core
- Thin materials, available in long-sole and short-sole versions
- Thin at the heel and edges, some available with a stabilizing heel cup

INFORMATION

MODIFICATIONS

- The core is partially or completely thermo-formable (hot air blower or oven, 100 to 120 °C, 30 to 60 sec. for partial thermo-forming, approx. 2 to 3 min. for complete thermo-forming)
- Easy to grind (grinding belt, grit 60–100, polishing brush on the blade disk)
- Can be glued to other materials (standard workshop adhesives)

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge (except leather top covers)
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

Short



TOP COVER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
–	–	Non-woven	Narrow, short, core visible	35–47	3 19541 1 00000**
–	–	Non-woven	Wide, short, core visible	35–48 ¹	3 19531 1 00000**
–	–	Non-woven	Wide, short, covered core	24–34 ¹⁺² 35–48	3 19231 1 00000**

With thermal reinforcement in the heel area



TOP COVER	METATARSAL PAD	BOTTOM COATING	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, semi-firm, 2 mm	Latex pad	–	Wide, short, thermal reinforcement heel	20–34 ¹⁺² 35–48	3 19151 1 00000**
Sleek leather, semi-firm, 2 mm	Latex pad	–	Wide, long, thermal reinforcement heel	20–34 ¹⁺² 35–48	3 19631 5 00000**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

GloboTec® Plastic / leather (Active)



ITEM NUMBER
3 10795 5 00000**

USES

- Correction and / or arch support for planovalgus foot
- Exceptionally stable material; also suitable for people with higher body weights and more significant deformities
- Suitable for a wide range of children's shoes

FEATURES

- Synthetic orthotic core
- Long and short-sole versions
- Thin at the heel and edges, with a stabilizing heel cup
- Loose supination wedges included

INFORMATION

MODIFICATIONS

- The core is partially or completely thermo-formable (hot air blower or oven, 120 to 130 °C, approx. 30 sec. for partial thermo-forming, approx. 2 to 3 min. for complete thermo-forming)
- Easy to grind (grinding belt, grit 40–100, polishing brush on the blade disk)
- Can be glued to other materials (standard workshop adhesives)

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Microfiber blue (animal prints)	–	✓	Short, high heel cup, loose supination wedges	20–42 ²	3 10790 1 00000**
–	Lunasoft (four colors) 2 mm (40° Shore A)	✓	Short, high heel cup, loose supination wedges	20–42 ²	3 10795 1 00000**
–	Lunasoft (four colors) 2 mm (40° Shore A)	✓	Long, high heel cup, loose supination wedges	20–42 ²	3 10795 5 00000**

² Double sizing

GloboTec® Globoplast



ITEM NUMBER

3 07973 1 00000**

USES

- For supporting flattened foot arches and non-surgical treatment of pigeon toe (metatarsus varus)
- Exceptionally stable material; also suitable for people with higher body weights and more significant deformities

FEATURES

- Synthetic orthotic core
- With a stabilizing heel cup

INFORMATION

MODIFICATIONS

- The core is partially or completely thermo-formable (hot air blower or oven, 100 to 180 °C, approx. 60 sec. for partial thermo-forming, approx. 4 to 6 min. for complete thermo-forming)
- Easy to grind (grinding belt, grit 24–100, polishing brush on the blade disk)
- Can be glued to other materials (standard workshop adhesives)

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

TOP COVER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	Metatarsal pad from size 36	¾-length, high heel cup, big toe strap	20–34 ² 36–46 ²	3 07973 1 00000**

² Double sizing



TRlactive® ONE STEP. MORE OPTIONS.

The new foot orthosis system allows the provision of premium products for all modifications.

The right product for all modifications					
		SUPPORTING EFFECT Control	CUSHIONING EFFECT Comfort	SOFT CUSHIONING EFFECT Soft	CUP-SHAPED Shape
TRlactive Select Line Page 45		•	•	•	
TRlactive Comfort Line Page 48			•	•	
TRlactive Soft Page 48				•	
TRlactive Shape Page 45		•			•

TRlactive® Select and TRlactive® Shape



TRlactive® Select

TRlactive® Select CHARGE-BASED 3-IN-1 PREMIUM SOLUTION

- + Orthotic blanks: require little modification and are easy to grind
- + Available in three versions, thanks to individual modification: supporting effect, cushioning effect, and soft cushioning effect
- + Thin design: suitable for a wide range of shoe types
- + Available in different versions with a supporting core and vertical functional zones

TRlactive®
SelectTRlactive®
Select heelTRlactive®
Select relaxTRlactive®
Select front relaxTRlactive®
Select rigidusTRlactive®
Select OATRlactive®
Shape

DESCRIPTION	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
TRlactive Select	TRlactive Technology	✓	Medium width	35–48	3 13500 5 00000**
	TRlactive Technology	✓	Wide	35–48	3 13850 5 00000**
TRlactive Select heel	TRlactive Technology	✓	Heel cushion	35–48	3 11500 5 00000**
TRlactive Select relax	TRlactive Technology	✓	Cushion long	35–48	3 136005 00000**
TRlactive Select front relax	TRlactive Technology	✓	Forefoot cushion	35–48	3 13700 5 00000**
TRlactive Select rigidus	TRlactive Technology	✓	Reinforced on the left	35–48	3 12600 5 00000**
	TRlactive Technology	✓	Reinforced on the right	35–48	3 12700 5 00000**
	TRlactive Technology	✓	Reinforced on both sides	35–48	3 12500 5 00000**
TRlactive Select OA	TRlactive Technology	✓	Supination	35–48	3 15600 5 00000**
	TRlactive Technology	✓	Pronation	35–48	3 14600 5 00000**
TRlactive Shape	TRlactive Technology	✓	Shell	35–48	3 10500 5 00000**

SOFT FOAM

SOFT FOAM ORTHOTIC BLANKS – FOR PRESSURE- SENSITIVE FEET



TRlactive® Comfort and TRlactive® Soft



- + Cushioning orthotic blanks in a comfort shoe shape: easy to grind
- + Vertical construction with three different materials
- + Can be adapted to suit a wide range of foot proportions and shoe shapes
- + Available in various versions



**TRlactive®
Comfort heel**



**TRlactive®
Comfort relax**



**TRlactive®
Comfort
front relax**



**TRlactive®
Comfort rigidus**



**TRlactive®
Soft**

DESCRIPTION	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
TRlactive Comfort heel	TRlactive Technology	✓	Heel cushion	35–50	3 24600 5 00000**
TRlactive Comfort relax	TRlactive Technology	✓	Cushion long	35–50	3 24100 5 00000**
TRlactive Comfort front relax	TRlactive Technology	✓	Forefoot cushion	35–50	3 24300 5 80000**
TRlactive Comfort rigidus	TRlactive Technology	✓	Reinforced on both sides	35–50	3 24700 5 00000**
TRlactive Soft	TRlactive Technology	✓	–	35–50	3 24250 5 00000**

ErgoPad® Soft style, complete+



ITEM NUMBER
3 82504 5 00008**

USES

- Orthotic blank for functional cushioning and relief
- Two versions for the individual treatment of foot problems
- For leisure and everyday shoes

FEATURES

- Soft cushion layer with absorbing functional zones (in the heel / midfoot and forefoot area, depending on the version)
- Lightweight product
- Can be adapted to suit a wide range of foot proportions and shoe shapes

INFORMATION

MODIFICATIONS

- Can be ground (lateral foot orthosis shape and foot orthosis thickness)
 - Easy to grind (grinding belt, grit 40–100, polishing brush on the blade disk)
 - Easy to glue (standard workshop adhesives)
 - Partly thermo-formable (hot air blower, up to 120 °C, partially 20–30 sec.)
- Please note: Material shrinkage possible!

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge (except leather top covers)
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

DESCRIPTION	TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
complete+	–	35 / 35° Shore A, gray	✓	Cushioned heel, base of 5th MTH, forefoot	35–48	3 82500 5 00008**
complete+	–	30 / 35° Shore A, orange	✓	Cushioned heel, base of 5th MTH, forefoot	35–48	3 82504 5 00008**

ErgoPad® system:s



ITEM NUMBER
3 28050 5 35864**

USES

- ErgoPad system:s offers a special base layer for every common foot misalignment:
- Base layer A: Mild arch decreased / splay foot
- Base layer B: Mild arch decreased foot
- Base layer C: Pronounced (valgus) arch decreased / splay foot
- Base layer D: Pronounced (valgus) arch decreased / splay foot with chronic heel and ankle pain
- Base layer E: Mild to moderate pes cavus / high arched foot

FEATURES

- Available in both unglued and glued versions
- Different top layer groups
- Standard shoe shape

INFORMATION

MODIFICATIONS

- Requires little modification
 - Easy to grind (grinding belt, grit 40–100, polishing brush on the blade disk)
 - Easy to glue (standard workshop adhesives)
 - Thermo-forming possible at approx. 80–100 °C, for approx. 30 sec.
- Please note: Material shrinkage possible!


CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge and / or a suitable disinfectant
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

THE DESIGN IN DETAIL

TOP LAYER

Selection of EVA / PE cushion layers with / without top cover



BASE LAYER



Base layers

DESCRIPTION	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Type A	Flat, anatomical	Long, thermal reinforcement, mild arch decreased / splay foot	35–50	3 28050 5 05640**
Type B	–	Long, thermal reinforcement, mild arch decreased foot	35–50	3 28050 5 05650**
Type C	Anatomical	Long, thermal reinforcement, pronounced (valgus) arch decreased / splay foot	35–50	3 28050 5 05660**
Type D	Anatomical	Long, thermal reinforcement, pronounced (valgus) arch decreased / splay foot with chronic heel and ankle pain	35–50	3 28050 5 05670**
Type E	Anatomical	Long, thermal reinforcement, mild to moderate pes cavus / high arched foot	35–50	3 28050 5 05680**

Glued versions, long, with thermal reinforcement

DESCRIPTION	TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Type A	Microfiber gray	TECfom light azure blue 3 mm (25° Shore A)	Flat, anatomical	Thermal reinforcement	35–50	3 28050 5 35806**
Type B	–	POWERDOTS pistachio 3 m	Sensorimotor spots	Thermal reinforcement	35–50	3 28050 5 35850**
Type C	–	TECfom light azure blue 3 mm (25° Shore A)	Anatomical	Thermal reinforcement	35–50	3 28050 5 35860**
Type C	–	TECtop silver-gray perforated 3 mm (30° Shore A)	Anatomical	Thermal reinforcement	35–50	3 28050 5 35864**
Type C	Microfiber gray	TECfom light azure blue 3 mm (25° Shore A)	Anatomical	Thermal reinforcement	35–50	3 28050 5 35866**
Type C	Microfiber brown	TECfom light azure blue 3 mm (25° Shore A)	Anatomical	Thermal reinforcement	35–50	3 28050 5 35868**
Type C	Microfiber gray	Globofit blue 3 mm (12° Shore A)	Anatomical	Thermal reinforcement	35–50	3 28050 5 35878**

Top layers

TOP COVER	CUSHION LAYER	METATARSAL PAD	SIZE (**)	ITEM NUMBER
Microfiber gray	TECfom, light azure blue 3 mm (25° Shore A)	–	35–50	3 28905 5 22051**
–	PowerDots pistachio 3 mm	Sensorimotor spots	35–50	3 68905 5 02562**

GloboTec® Soft




GloboTec® Soft
microfiber soft-cushioning
orthotic blank

ITEM NUMBER
3 22770 5 01400**

THE DESIGN IN DETAIL

- 1 TOP COVER
- 2 CUSHION LAYER
- 3 METATARSAL PAD
- 4 THERMOPLASTIC REINFORCEMENT
- 5 BASE LAYER



USES

- Soft foam orthotic blank for functional cushioning and relief
- For treatment of various foot deformities
- Depending on the cushion layer / base layer combination, suitable e.g. for sport and leisure activities, treating pressure-sensitive feet and in the case of diabetic foot syndrome or rheumatoid foot deformities
- For leisure and comfort shoes

FEATURES

- Adequate orthopedic support for the foot arches
- Can be adapted to suit a wide range of foot proportions and shoe shapes
- Base layer and cushion layers available with different levels of density

INFORMATION

MODIFICATIONS

- Must be ground (foot orthosis shape at the sides and foot orthosis thickness)
 - Easy to grind (grinding belt, grit 40–100, polishing brush on the blade disk)
 - Easy to glue (standard workshop adhesives)
 - Partially or completely thermo-formable (hot air blower or oven, up to 120 °C, approx. 20–30 sec. for partial thermo-forming, approx. 2 min. for complete thermo-forming)
- Please note: Material shrinkage possible!

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge (except leather top covers)
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

Soft foam blanks, single-layer TECfom (25° Shore A), long



TOP COVER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	✓	Narrow	34 – 48	3 22030 5 00000**
–	✓	Wide	34 – 48	3 22050 5 00000**
Calfskin partly perforated	✓	Wide	35 – 48	3 22130 5 00000**
Calfskin partly perforated	✓	Wide	34 – 48	3 22150 5 00000**
Microfiber beige	✓	Wide	34 – 48	3 22150 5 01400**

Soft foam blanks, single-layer TECmoll (35° Shore A), long



TOP COVER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	✓	Comfort shoe shape	35 – 48	3 22200 5 00000**
–	✓	Wide	20 – 34 ¹⁺² 35 – 48	3 22250 5 00000**
–	–	Wide	35 – 48	3 22250 5 00350**
–	✓ from size 35	Wide, heel cup	20 – 34 ² 35 – 48	3 22270 5 00000**

Supporting soft foam blanks, single-layer TECstrong (50° Shore A)



TOP COVER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
onSteam crema partly perforated	✓	Thermal reinforcement	35 – 48	3 21255 5 00000**
–	✓	Thermal reinforcement	35 – 48	3 21260 5 00000**

Soft foam blanks, TECsoft / TECsoft (25° / 25° Shore A), long



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECsoft gray 5 mm	Latex pad	Medium width, thermal reinforcement	35 – 48	3 20230 5 00000**
–	TECsoft gray 5 mm	Latex pad	Wide, thermal reinforcement	26 – 34 ² 35 – 50	3 20250 5 00000**
Calfskin partly perforated	TECsoft gray 5 mm	Latex pad	Medium width, thermal reinforcement	35 – 48	3 20330 5 00000**
Calfskin partly perforated	TECsoft gray 5 mm	Latex pad	Wide, thermal reinforcement	24 – 34 ¹⁺² 35 – 50	3 20350 5 00000**
Calfskin partly perforated	TECsoft gray 3 mm	Latex pad	Wide, thermal reinforcement	26 – 34 ¹⁺² 35 – 49	3 20530 5 00000**
–	TECsoft gray 3 mm	✓	Wide, thermal reinforcement	35 – 42	3 24702 5 00000**
–	TECsoft gray 3 mm	✓	Wide, thermal reinforcement, square shape	38 – 48	3 24703 5 00000**

Soft foam blanks, TECsoft / TECmoll (25° / 35° Shore A), long



TOP COVER	CUSHION LAYER	BASE LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECsoft black 5 mm	TECmoll 35° Shore A	✓	Wide, thermal reinforcement	35 – 48	3 20250 5 10000**
–	TECsoft black 5 mm	TECmoll 35° Shore A	✓	Wide, thermal reinforcement, comfort shoe shape	35 – 48	3 20230 5 10100**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

GloboTec® Soft



Soft foam blanks, TECfom / TECmoll (25° / 35° Shore A), long

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECfom gray 6 mm	✓	Wide	20–34 ² 35–50	3 20550 5 00000**
–	TECfom gray 6 mm	✓	Wide, thermal reinforcement	20–34 ² 35–50	3 20550 5 04721**
Microfiber beige	TECfom gray 6 mm	✓	Wide	20–34 ² 35–50	3 20650 5 01400**
onSteam anthracite partly perforated	TECfom gray 6 mm	✓	Wide	35–48	3 20650 5 01422**
Microfiber beige	TECfom gray 6 mm	✓	Wide, thermal reinforcement	20–34 ¹⁺² 35–50 ¹	3 20650 5 99772**
–	TECfom blue 6 mm	✓	Wide	20–34 ² 35–48	3 20850 5 00000**
Calfskin partly perforated	TECfom blue 6 mm	✓	Wide	20–34 ² 35–48	3 20950 5 00000**
Microfiber blue	TECfom blue 6 mm	✓	Wide	20–34 ¹⁺² 35–48	3 20950 5 01401**



Soft foam blanks, Globoflex / TEClast (15° / 45° Shore A), long

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	Globoflex 6 mm, PUR top, PE bottom	✓	Thermal reinforcement, comfort shoe shape	35–48	3 22700 5 00000**
–	Globoflex 6 mm, PE top, PUR bottom	✓	Thermal reinforcement, comfort shoe shape	35–48	3 22710 5 00000**
Microfiber beige	Globoflex 6 mm, PUR top, PE bottom	✓	Thermal reinforcement, comfort shoe shape	35–48	3 22720 5 01400**
–	Globoflex 6 mm, PUR top, PE bottom	✓	Wide, thermal reinforcement	20–34 ¹⁺² 35–48	3 22750 5 00000**
Microfiber beige	Globoflex 6 mm, PUR top, PE bottom	✓	Wide, thermal reinforcement	20–34 ² 35–50	3 22770 5 01400**
Microfiber blue	Globoflex 6 mm, PUR top, PE bottom	✓	Wide, thermal reinforcement	20–34 ¹⁺² 35–50	3 22770 5 01401**



Soft foam blanks, Globofit / TEClast (12° / 45° Shore A), long

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Calfskin partly perforated	Globofit beige 3 mm	Latex pad	Wide	20–34 ¹⁺² 35–48	3 22650 5 00000**
Calfskin partly perforated	Globofit beige 3 mm	Latex pad	Wide, thermal reinforcement	20–34 ¹⁺² 35–48	3 22650 5 04722**



Soft foam blanks, Poron Kombi / TEClast (– / 45° Shore A), long

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	Poron Kombi 6 mm	✓	Comfort shoe shape	35–48	3 21500 5 02296**
–	Poron Kombi 6 mm	✓	Comfort shoe shape, thermal reinforcement	35–48	3 21500 5 12296**
–	Poron Kombi 6 mm	✓	Wide	26–34 ¹⁺² 35–50	3 21550 5 02296**
–	Poron Kombi 6 mm	✓	Wide, thermal reinforcement	26–34 ¹⁺² 35–50	3 21550 5 12296**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

Soft foam blanks, Poron Kombi / TECmoll / TEClast (15° / 35° / 45° Shore A), long


TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	Poron Kombi 7.5 mm	–	Comfort shoe shape, thermal reinforcement	35–50	3 23800 5 00000**
–	Poron Kombi 7.5 mm	–	Comfort shoe shape, thermal reinforcement, loose cover	35–50	3 23800 5 00440**

Soft foam blanks, TECtop / TEClast (30° / 45° Shore A), long


TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECtop blue perforated 4 mm	✓	Wide	20–34 ² 35–50	3 21550 5 02122**
–	TECtop brown perforated 4 mm	✓	Wide	35–48	3 21550 5 02152**
–	TECtop blue perforated 4 mm	✓	Wide, thermal reinforcement	20–34 ¹⁺² 35–50	3 21550 5 11266**
–	TECtop blue perforated 4 mm	✓	Wide	35–48	3 22850 5 02122**

Soft foam blanks, TECmoll / TEClast (35° / 45° Shore A), long


TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECmoll blue 6 mm	✓	Comfort shoe shape	35–48	3 21800 5 00000**
–	TECmoll blue 6 mm	✓	Wide	20–34 ² 35–50	3 21850 5 00000**
–	TECmoll blue 6 mm	–	Wide	20–34 ² 35–48	3 21850 5 00350**
Microfiber blue	TECmoll blue 6 mm	✓	Wide	20–34 ² 35–50	3 21950 5 01401**
–	TECmoll gray 6 mm	✓	Comfort shoe shape	35–48	3 21500 5 00000**
–	TECmoll gray 6 mm	✓	Wide	20–34 ² 35–50	3 21550 5 00000**
–	TECmoll gray 6 mm	–	Wide	35–50	3 21550 5 00350**
–	TECmoll gray 6 mm	✓	Wide, thermal reinforcement	35–50	3 21550 5 04722**
–	TECmoll gray 6 mm	✓ From size 35	Wide, heel cup	20–34 ² 35–48	3 21570 5 00000**
Calfskin partly perforated	TECmoll gray 6 mm	✓	Comfort shoe shape	35–48	3 21600 5 00000**
Microfiber beige	TECmoll gray 6 mm	✓	Comfort shoe shape	35–48	3 21600 5 01400**
Calfskin partly perforated	TECmoll gray 6 mm	✓	Wide	20–34 ² 35–50	3 21650 5 00000**
Microfiber beige	TECmoll gray 6 mm	✓	Wide	20–34 35–48	3 21650 5 01400**
onSteam anthracite partly perforated	TECmoll gray 6 mm	✓	Wide	35–48	3 21650 5 01422**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

GloboTec® Soft

Soft foam blanks, TECtop / thermal cork (30° / 50° Shore A), long



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECtop blue perforated 4 mm	✓	Narrow, thermal reinforcement	24–34 ² 34–50	3 58650 5 04722**
–	TECtop blue perforated 4 mm	✓	Wide	24–32 ² 34–48	3 58650 5 00000**
–	TECtop blue perforated 4 mm	✓ From size 36	Wide, thermal reinforcement, high heel cup	20–34 ² 36–48 ²	3 58890 5 00000**

Soft foam blanks, TECtop / TECstrong (30° / 55° Shore A), long



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECtop blue 4 mm	✓	Narrow, thermal reinforcement, square shape	35–42	3 50552 5 80000**
–	TECtop blue 4 mm	✓	Wide, thermal reinforcement, square shape	35–48	3 50553 5 80000**
–	TECtop blue perforated 4 mm	✓	Wide, thermal reinforcement	35–42	3 24662 5 02122**
–	TECtop red perforated 4 mm	✓	Wide, thermal reinforcement	35–42	3 24662 5 02132**
–	TECtop blue perforated 4 mm	✓	Wide, thermal reinforcement, square shape	38–48	3 24663 5 02122**

Supporting soft foam blanks, TECstrong / TECstrong (55° / 55° Shore A), long



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECstrong brown 6 mm	✓	Wide	35–42	3 24772 5 00000**
–	TECstrong brown 6 mm	✓	Wide, square shape	38–48	3 24773 5 00000**

Soft foam blanks with detorsion cut, TECtop / TECtop (30° / 30° Shore A), long



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECtop white 6 mm	✓	Wide, thermal reinforcement	35–42	3 24742 5 00000**
–	TECtop white 6 mm	✓	Wide, thermal reinforcement, square shape	36–48	3 24743 5 00000**

² Double sizing



Soft foam blanks with detorsion cut, TECtop / TEClast (30° / 45° Shore A), long

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECtop white 6 mm	✓	Wide	35–42	3 24762 5 00000**
–	TECtop white 6 mm	✓	Wide, square shape	36–48	3 24763 5 00000**



Soft foam blanks with heel cup, long

TOP COVER	CUSHION LAYER	BASE LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TEClast gray 6 mm (45° Shore A)	TEClast gray 6 mm (45° Shore A)	✓ From size 35	Wide, heel cup	20–34 ² 35–48	3 21270 5 00000**
–	TECtop blue perforated 4 mm (30° Shore A)	TEClast blue 6 mm (45° Shore A)	✓ From size 35	Wide, high heel cup	20–34 ² 36–46 ²	3 21890 5 00000**
–	Lunasoft (four colors) 3 mm (40° Shore A)	TECstrong blue 10 mm (55° Shore A)	–	Wide, high heel cup	20–34 ² 36–40 ²	3 24681 5 00000**



Other soft foam blanks, long

TOP COVER	CUSHION LAYER	BASE LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECtop white 6 mm (30° Shore A)	TECstrong white 10 mm (55° Shore A)	✓	Poron heel cushion, comfort shoe shape	35–48	3 20400 5 84152**
–	TECtop white 6 mm (30° Shore A)	TECstrong white 10 mm (55° Shore A)	✓	Poron heel cushion, thermal reinforcement, comfort shoe shape	35–48	3 20400 5 85152**
–	TEClast gray 6 mm (45° Shore A)	TEClast gray 6 mm (45° Shore A)	✓	Wide	20–34 ¹⁺² 35–48	3 21250 5 00000**
–	TECfom gray 6 mm (25° Shore A)	TECstrong white 6 mm (45° Shore A)	✓	Comfort shoe shape, thermal reinforcement	35–40	3 24612 5 01999**
–	Bitec (25°–35° Shore)	TECstrong white 10 mm (45° Shore A)	✓	Wide, standard heel cup, thermal reinforcement	35–48	3 59000 5 80000**



Berkosoft Vario base layers

TOP COVER	CUSHION LAYER	METATARSAL PAD	BASE LAYER	VERSION	SIZE (**)	ITEM NUMBER
–	–	–	EVA soft foam blue	Wide, short, semi-firm	36–40 ²	3 24912 1 00000**
–	–	–	EVA soft foam blue	Wide, short, semi-firm	42–48 ²	3 24913 1 00000**
–	–	–	EVA soft foam green	Wide, short, soft	36–40 ²	3 24922 1 00000**
–	–	–	EVA soft foam green	Wide, short, soft	42–48 ²	3 24923 1 00000**
–	–	–	Thermal cork	Wide, short, firm	36–40 ²	3 24932 1 00000**
–	–	–	Thermal cork	Wide, short, firm	42–48 ²	3 24933 1 00000**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

GloboTec® Soft



Berkosoft top layers

TOP COVER	CUSHION LAYER	METATARSAL PAD	BASE LAYER	VERSION	SIZE (**)	ITEM NUMBER
Colored fabric Cambrelle	PE cushion layer	–	–	Wide, long, very soft	36–40 ²	3 24972 5 00000**
Colored fabric Cambrelle	PE cushion layer	–	–	Wide, long, very soft	42–48 ²	3 24973 5 00000**
–	EVA cushion layer yellow	–	–	Wide, long, soft	36–40 ²	3 24982 5 00000**
Microfiber	EVA cushion layer yellow	–	–	Wide, long, soft	36–40 ²	3 24982 5 01220**
–	EVA cushion layer yellow	–	–	Wide, long, soft	42–48 ²	3 24983 5 00000**
Microfiber	EVA cushion layer yellow	–	–	Wide, long, soft	42–48 ²	3 24983 5 01220**
–	EVA cushion layer blue	–	–	Wide, long, soft	36–40 ²	3 24992 5 00000**
Microfiber	EVA cushion layer blue	–	–	Wide, long, soft	36–40 ²	3 24992 5 01220**
–	EVA cushion layer blue	–	–	Wide, long, soft	42–48 ²	3 24993 5 00000**
Microfiber	EVA cushion layer blue	–	–	Wide, long, soft	42–48 ²	3 24993 5 01220**

² Double sizing



Soft foam blanks for children, long

TOP COVER	CUSHION LAYER	BASE LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
–	TECsoft gray 5 mm (25° Shore A)	TECsoft gray 5 mm (25° Shore A)	Latex pad	Thermal reinforcement	24–34 ²	3 20250 5 00000**
Calfskin partly perforated	TECsoft gray 5 mm (25° Shore A)	TECsoft gray 5 mm (25° Shore A)	Latex pad	Thermal reinforcement	24–34 ¹⁺²	3 20350 5 00000**
Calfskin partly perforated	TECsoft gray 5 mm (25° Shore A)	TECsoft gray 5 mm (25° Shore A)	Latex pad	Thermal reinforcement	26–34 ¹⁺²	3 20530 5 00000**
–	Lunasoft (four colors) 3 mm (40° Shore A)	TECtop white 6 mm (30° Shore A)	–		20–34 ²	3 24741 5 02241**
–	TECfom gray 6 mm (25° Shore A)	TECmoll gray 6 mm (35° Shore A)	✓		20–34 ²	3 20550 5 00000**
Microfiber beige	TECfom gray 6 mm (25° Shore A)	TECmoll gray 6 mm (35° Shore A)	✓		20–34 ²	3 20650 5 01400**
–	TECfom blue 6 mm (25° Shore A)	TECmoll blue 6 mm (35° Shore A)	✓		20–34 ²	3 20850 5 00000**
–	Globoflex PUR top, PE bottom 6 mm (15° Shore A)	TEClast gray 6 mm (45° Shore A)	✓	Thermal reinforcement	20–34 ¹⁺²	3 22750 5 00000**
Microfiber beige	Globoflex PUR top, PE bottom 6 mm (15° Shore A)	TEClast gray 6 mm (45° Shore A)	✓	Thermal reinforcement	20–34 ²	3 22770 5 01400**
Microfiber blue	Globoflex PUR top, PE bottom 6 mm (15° Shore A)	TEClast gray 6 mm (45° Shore A)	✓	Thermal reinforcement	20–34 ¹⁺²	3 22770 5 01401**
Calfskin partly perforated	Globofit beige 3 mm (12° Shore A)	TEClast gray 6 mm (45° Shore A)	Latex pad		20–34 ¹⁺²	3 22650 5 00000**
Calfskin partly perforated	Globofit beige 3 mm (12° Shore A)	TEClast gray 6 mm (45° Shore A)	Latex pad	Thermal reinforcement	20–34 ¹⁺²	3 22650 5 04722**
–	Poron Kombi 6 mm	TEClast gray 6 mm (45° Shore A)	✓		26–34 ¹⁺²	3 21550 5 02296**
–	TECtop blue perforated 4 mm (30° Shore A)	Thermal cork 8 mm (50° Shore A)	✓ From size 36	Thermal reinforcement, high heel cup	20–34 ²	3 58890 5 00000**
–	TEClast gray 6 mm (45° Shore A)	TEClast gray 6 mm (45° Shore A)	✓ From size 35	Heel cup	20–34 ²	3 21270 5 00000**
–	TECtop blue perforated 4 mm (30° Shore A)	TEClast blue 6 mm (45° Shore A)	✓ From size 35	High heel cup	20–34 ²	3 21890 5 00000**
–	Lunasoft (four colors) 3 mm (40° Shore A)	TECstrong blue 10 mm (55° Shore A)	–	High heel cup	20–34 ²	3 24681 5 00000**
–	TEClast gray 6 mm (45° Shore A)	TEClast gray 6 mm (45° Shore A)	✓		20–34 ¹⁺²	3 21250 5 00000**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

MILLED ORTHOSES

MILLED ORTHOTIC BLANKS

PERFECT
CUSTOM-MADE
FOOT PRODUCTS

BODYTRONIC® ID:CAM

YOUR PERSONAL TOOLBOX

Turn individual customer requests into reality with Bodytronic ID:CAM.

Using the BODYTRONIC ID:CAM system solution, you can model customized milled orthoses that offer a precise fit.

SELECT THE PROFILE THAT SUITS YOU



BASIC BAUERFEIND FULL SERVICE

- You send the necessary customer data and information concerning the material, shape and shoe size.
- Bauerfeind takes over the modeling and manufacture of the milled orthoses for you.





ADVANCED BAUERFEIND VALIDATION SERVICE

- You carry out the modeling process based on the customer data using an integrated model and module library.
- Bauerfeind manufactures your milled orthoses after verification.



PROFESSIONAL BAUERFEIND MANUFACTURING SERVICE

- You perform the configuration and modeling yourself and then send the order online directly into the production.
- Bauerfeind manufactures your custom-milled orthoses without checking them beforehand.




The free web application is available 24 / 7 at

idcamrev2.bauerfeind.com




YOUR BENEFITS AT A GLANCE

- **PERSONALIZED CONSULTATION ADAPTED TO PATIENT NEEDS**
using comprehensive digital medical history form
- **INDIVIDUALIZED MODELING PROCESS**
thanks to user-based template library
- **NOTICEABLE TIME SAVINGS**
thanks to clear and intuitive user interface as well as personal library
- **FAST AND TAILORED CUSTOMER CARE**
thanks to foot orthosis delivery within 48 hours
- **COMPREHENSIVE OPTIONS FOR HIGHLY SPECIALIZED PRODUCTS**
thanks to different modeling modules for the treatment of adults, children and patient with diabetes, including plus sizes and customized milled foot orthoses with core element
- **MODELING OF DIABETES-ADAPTED FOOT CUSHIONING by using**
a digital last in the diabetes module
- **USE OF ALL CONVENTIONAL RECORDING METHODS**
based on the upload of different data (2D, 3D, foot pressure measurement)
- **MATERIAL SELECTION ADAPTED TO INDIVIDUAL CUSTOMER NEEDS**
thanks to a huge choice of single-, dual- and triple-layer materials with different degrees of density



MEASURE – MODEL – VALIDATE

You can decide which stages of the foot orthosis production process you want to participate in. If you like, Bauerfeind can take the digital measurement data and start modeling and manufacturing the milled orthoses straight away.



ErgoPad® ID:



ErgoPad® ID:
Trend dynamic

USES

- BODYTRONIC ID:CAM is Bauerfeind's concept for manufacturing precisely milled customized foot orthoses to suit your specific indication and foot care requirements.
- Orthoses are created using an intuitive and free web application offering three different user profiles
- Based on 2D data (scan, blueprint), foot pressure measurement data and 3D data (imprint foam scan)
- Suitable for comfort, leisure and sports shoes, depending on the configuration of the foot orthosis

FEATURES

- Thanks to the wide variety of material combinations available, a solution can be found for every foot problem
- The integrated and freely modifiable model and module library helps you create customized milled foot orthoses quickly and efficiently
- Manufactured using precise CNC milling machines and deliverable / delivered within a short time
- Forefoot and / or heel intarsia (cushions) can be incorporated into the foot orthosis
- Foot orthoses are also available with a top cover (microfiber / lambskin leather)
- Further processing by an orthotist is required
- Can be delivered pre-ground at the sides

INFORMATION

MODIFICATIONS

- Must be ground (to fit the shoe / side of the foot if required)
- Easy to grind (grinding belt, grit 60–100, polishing brush on the blade disk)
- Easy to glue (standard workshop adhesives)
- Partially thermo-formable (hot air blower, up to 120 °C, 20–30 sec.)
- Avoid attempting large-scale thermoforming
- A top cover, e.g. microfiber, is recommended for mechanical protection of the soft foam material

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge (except leather top covers)
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

DESCRIPTION		INNER SOLE SHAPE	BASE THICKNESS	FOREFOOT MODULES	MIDFOOT MODULES	SIZE
ID:Trend	With its moderate arch supports , this predesign is a perfect starting point for treating high-arched / arch decreased feet or splay feet	For comfort shoes	3.5 mm	None	Metatarsal pad (long), medial arch support (long)	35–50
ID:Trend Dynamic	With its proprioceptive modules , this predesign provides a good basis for stimulating the foot muscles and thus actively correcting misalignments	For comfort shoes	3 mm	Toe bank	Metatarsal pad (large), medial arch support (long), lateral arch support (short)	35–50
ID:Style	The moderate arch supports are the perfect starting point for treating high-arched / arch decreased feet or splay feet; preselected lateral arch supports allow easy adaption to combat excessive supination; the thin forefoot area makes shortening of the foot orthosis easy	For comfort shoes	1.5 mm	None	Metatarsal pad (long), medial arch support (short), lateral arch supports (long)	35–50
ID:Style Dynamic	With its proprioceptive modules , this predesign provides a good basis for stimulating the foot muscles and thus actively correcting misalignments; thanks to its low forefoot height , it's perfect for shoes with little space in the toe area	For comfort shoes	1.5 mm	Toe bank	Metatarsal pad (large), medial arch support (short), lateral arch support (short)	35–50
ID:Fit	The moderate arch supports are a perfect starting point for treating high-arched / arch decreased feet or splay feet	For standard shoes	3.5 mm	None	Metatarsal pad (long), medial arch support (long)	35–50
ID:Fan Cup	The high heel cup holds the heel in a secure position and is a good basis for treating high-arched / arch decreased feet or splay feet	For standard shoes	2.5 mm	None	Metatarsal pad (long), heel cup (long)	35–50
ID:Bodytronic 1	Slim feet with problems relating to high-arched, arch decreased or splay feet can be treated particularly well using moderate arch supports	Standard shoe shape, narrow	1.5 mm	None	Metatarsal pad (long), medial arch support (long)	35–48
ID:Bodytronic 2	Average-width feet with problems relating to high-arched, arch decreased or splay feet can be treated well using moderate arch supports	Standard shoe shape, medium width	1.5 mm	None	Metatarsal pad (long), medial arch support (long)	35–48
ID:Bodytronic 3	Wide feet with problems relating to high-arched, arch decreased or splay feet can be treated in a targeted manner using moderate arch supports	Comfort shoe shape, wide	1.5 mm	None	Metatarsal pad (long), medial arch support (long)	35–48
ID:Extra-wide	This predesign has been specially developed for extra-wide and high-volume feet ; moderate valgus and arch decreased foot deformities can be treated in a targeted way; additionally, the wide contour allows fitting them into therapeutic shoes	Extra wide	1.5 mm	None	Metatarsal pad (long), medial arch support (long)	35–48

ErgoPad® ID:

DESCRIPTION		INNER SOLE SHAPE	BASE THICKNESS	FOREFOOT MODULES	MIDFOOT MODULES	SIZE
ID:Trend Flat	With this predesign, this foot orthosis can be modeled completely without a predefined topology and thus tailored to customer needs	For comfort shoes	1.5 mm	None	None	35 – 50
ID:Toy	A predesign for the treatment of children with a high, extended heel cup and toe grip; a predesign is also available for treating children – ErgoPad ID: Toy Flat – which makes it possible to design a foot orthosis entirely without any predefined topology	For children's shoes	2 mm	Toe grip (long)	Metatarsal pad (long), heel cup (extended, high)	21 – 34
Diabetes 1	With its moderate arch support , this predesign is a perfect starting point for treating high-arched / arch decreased feet or splay feet	Comfort shoe shape (freely modifiable)	3.5 mm	None	Metatarsal pad (long), medial arch support (long)	35 – 50
Diabetes 2	The moderate arch supports are the perfect starting point for treating high-arched / arch decreased feet or splay feet; preselected lateral arch supports allow easy adaption to combat excessive supination; the thin forefoot area makes shortening of the foot orthosis easy	Comfort shoe shape (freely modifiable)	1.5 mm	None	Metatarsal pad (long), medial arch support (short), lateral arch support (long)	35 – 50
Diabetes 3	With its moderate modules , this predesign is a perfect starting point for treating high-arched / arch decreased feet or splay feet	Comfort shoe shape (freely modifiable)	3.5 mm	None	Metatarsal pad (long), medial arch support (long)	35 – 50
Diabetes 4	Thanks to its moderate modules , this pre-design provides a good basis for treating sensitive and injury-prone feet ; a pre-selected moderate heel cup also ensures good foot guidance	Comfort shoe shape (freely modifiable)	8 mm	None	Medial arch support (short), heel cup (long)	35 – 50

Additional features

DESCRIPTION	VERSION
Top cover	onSteam / microfiber / lambskin leather
Heel cushion	9.5 mm cushioning material
Forefoot / hindfoot intarsia	6 mm cushion material

MILLED SPORTS FOOT ORTHOSES



FOOTBALL

ITEM NUMBER
3 75500 0 00000**



GOLF RIGHT-HANDED SWING

ITEM NUMBER
3 75100 0 00000**

USES

- Designed to meet sport-specific requirements
- Biomechanically optimized support, stability and cushioning for the feet, knees, hip joints and spine
- Suitable for golf, inline skating, mountain biking, Alpine skiing and football
- Suitable for sports shoes

FEATURES

- CNC-milled, multi-layered, horizontal soft foam blank
- Sport-specific design in accordance with biomechanical requirements
- Orthopedically effective support for the foot arches
- Can be adapted to suit a wide range of foot proportions and shoe shapes

INFORMATION

MODIFICATIONS

- Must be ground (to fit the shoe / side of the foot)
- Easy to grind (grinding belt, grit 60–100, polishing brush on the blade disk)
- Easy to glue (standard workshop adhesives)
- Top cover recommended, e.g. high-quality microfiber

CLEANING / CARE

- Surfaces can be wiped clean using a mild cleanser and a sponge
- Not machine-washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

DESCRIPTION	VERSION	SIZE (**)	ITEM NUMBER
Golf right-handed swing	35 / 25 / 50° Shore A	35–48	3 75100 0 00000**
Golf left-handed swing	35 / 25 / 50° Shore A	35–48 ¹	3 75150 0 00000**
Skating	35 / 25 / 50° Shore A	35–48 ¹	3 75200 0 00000**
Alpine skiing	35 / 25 / 50° Shore A	35–48	3 75300 0 00000**
Mountain biking	35 / 25 / 50° Shore A	35–48 ¹	3 75400 0 00000**
Football	45 / 25° Shore A + base material	35–48	3 75500 0 00000**

¹ Not a stock item. Please note that a longer delivery time will apply.

THE BEST OF TWO WORLDS – WEIGHTFLEX MEETS ID:CAM

The new generation of ErgoPad ID:CoreX orthoses combines a thin, self-supporting core with a customized, digitally modeled soft cushion surface.




ErgoPad® ID:CoreX

A milled foot orthosis with weightflex core

Precisely modeled
cushion layer with self-
supporting core

Individual modeling in the
BODYTRONIC ID:CAM web
application



Core available in soft, medium,
and strong for different load
requirements

Comprehensive
information relating to
weightflex technology
can be found on page 13.



Synthetic core with
weightflex technology

SIZE


35–50

Can be modeled and ordered via
ID:CAM web application.

More information:



ONLY BAREFOOT IS BETTER



Discover now! →

ErgoPad® ID:SANDAL

Enjoying the summer with
custom-made sandals

Summer sandal with an
orthopedic footbed

Individual modeling in the
BODYTRONIC ID:CAM web
application



Many colors for straps and
cushion layers

A one-off with the perfect fit,
thanks to CAD milling technology

STRAP COLORS



CUSHION LAYER COLORS



Delivered as an assembly kit

SIZE

35–50

Can be modeled and ordered via
ID:CAM web application.

More information:



CORK

CORK ORTHOTIC BLANKS FROM NATURE TO YOUR SHOES


GloboTec® Elastra-Cork

ITEM NUMBER
3 02522 5 00000**



THE DESIGN IN DETAIL

- 1 TOP COVER
- 2 CUSHION LAYER
- 3 METATARSAL PAD
- 4 THERMAL REINFORCEMENT
- 5 CORK BASE LAYER



USES

- Cork (and leather) orthotic blank to provide stable support for the foot arches
- Suitable for the treatment of typical foot deformities, e.g. arch decreased / splay foot
- A cushioning effect is also possible, depending on the type of cushion layer used
- Suitable for leisure and comfort shoes

FEATURES

- Bauerfeind's proprietary combination of natural rubber and cork is lightweight, elastic and dimensionally stable
- Can be adapted to suit a wide range of foot proportions and shoe shapes
- Long-sole and short-sole blanks available
- Cork base available in various thicknesses and shapes, with or without a wooden spline
- High-quality, vegetable tanned sleek leather ensures better conditions within the shoe

INFORMATION

MODIFICATIONS

- Must be ground (foot orthosis shape at the sides and foot orthosis thickness)
- Easy to grind (grinding belt, grit 40–100, polishing brush on the blade disk)
- Easy to glue (standard workshop adhesives)

CLEANING / CARE

- Not washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months



Base layers, short

METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
✓	Narrow, without reinforcement	35–42	3 02602 1 00000**
✓	Narrow, thin, without reinforcement	35–42	3 02502 1 00000**
✓	Narrow, wood veneer	35–42	3 02612 1 00000**
✓	Narrow, thick, wood veneer	35–42	3 02562 1 00000**
✓	Wide, without reinforcement	38–48	3 02603 1 00000**
✓	Wide, thick, without reinforcement	38–48	3 02623 1 00000**
✓	Wide, wood veneer	38–48	3 02613 1 00000**
✓	Wide, thick, wood veneer	38–48 ¹	3 02563 1 00000**

Without reinforcement, wide



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 2 mm	–	✓	Short	38–48	3 02163 1 00000**
Sleek leather, soft 2 mm	–	✓	Long, square shape	38–48	3 02163 5 00000**
–	TECsoft gray 5 mm (25° Shore A)	✓	Long, square shape	38–48	3 02523 5 00000**
Goatskin 0.8 mm	TECsoft gray 5 mm (25° Shore A)	✓	Long, square shape	38–48	3 02543 5 00000**

Without reinforcement, narrow



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 2 mm	–	✓	Short	20–24 ² 35–44	3 02162 1 00000**
Sleek leather, soft 2 mm	–	✓	Long	35–44	3 02162 5 00000**
–	TECsoft gray 5 mm (25° Shore A)	✓	Long	35–42	3 02522 5 00000**
Goatskin 0.8 mm	TECsoft gray 5 mm (25° Shore A)	✓	Long	35–42	3 02542 5 00000**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

GloboTec® Elastra-Cork

With wooden spline, wide



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 2 mm	–	✓	Long, thin, square shape	35–48	3 02023 5 04601**
Sleek leather, soft 2 mm	–	✓	Short	38–48	3 02163 1 04601**
Sleek leather, soft 2 mm	–	✓	Long, square shape	37–48	3 02163 5 04601**
Sleek leather, soft 2 mm	–	✓	Long, thick, square shape	40–50	3 02203 5 04601**
–	TECsoft gray 5 mm (25° Shore A)	✓	Long, square shape	38–50	3 02523 5 04601**
Goatskin 0.8 mm	TECsoft gray 5 mm (25° Shore A)	✓	Long, square shape	38–48	3 02543 5 04601**

With wooden spline, narrow



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 2 mm	–	✓	Long, thin	35–42	3 02022 5 04601**
Sleek leather, soft 2 mm	–	✓	Short	35–44	3 02162 1 04601**
Sleek leather, soft 2 mm	–	✓	Long	35–44	3 02162 5 04601**
–	TECsoft gray 5 mm (25° Shore A)	✓	Long	35–44	3 02522 5 04601**
Goatskin 0.8 mm	TECsoft gray 5 mm (25° Shore A)	✓	Long	35–42	3 02542 5 04601**

With thermal reinforcement, wide



TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 2 mm	–	✓	Long, square shape	38–48 ¹	3 02163 5 04721**
Sleek leather, soft 2 mm	–	Latex pad	Long, square shape	38–48	3 02163 5 22211**
Sleek leather, soft 2 mm	–	✓	Long, thick, square shape	35–50	3 02203 5 00000**
Goatskin 0.8 mm	TECsoft gray 5 mm (25° Shore A)	✓	Long, square shape	38–48	3 02543 5 04721**
Sleek leather, soft 1.25 mm	TECfom gray 3 mm (25° Shore A)	✓	Long	38–48 ¹	3 02543 5 84243**
Cowlining 1 mm	Globofit beige 2 mm (12° Shore A)	✓	Long	38–48 ¹	3 02543 5 84247**

¹ Not a stock item. Please note that a longer delivery time will apply.

With thermal reinforcement, narrow


TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 2 mm	–	✓	Long	35–44	3 02162 5 04721**
Sleek leather, soft 2 mm	–	Latex pad	Long	35–44	3 02162 5 22211**
Sleek leather, soft 2 mm	–	✓	Long, thick	35–43	3 02202 5 00000**
Goatskin 0.8 mm	TECsoft gray 5 mm (25° Shore A)	✓	Long	35–42 ¹	3 02542 5 04721**
Sleek leather, soft 1.25 mm	TECfom gray 3 mm (25° Shore A)	✓	Long	38–48 ¹	3 02542 5 84243**
Cowlining 1 mm	Globofit beige 2 mm (12° Shore A)	✓	Long	35–43 ¹	3 02542 5 84247**

Children's foot orthoses with a cupped design


TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 2 mm	–	✓	Long, without reinforcement	20–32 ² 33–34	3 02161 5 00000**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

GloboTec® Alpha-Cork



ITEM NUMBER
3 01650 5 04721**

USES

- Cork / leather orthotic blank to provide stable support for the foot arches
- Suitable for the treatment of typical foot deformities, e.g. arch decreased / splay foot
- A cushioning effect is also possible, depending on the type of cushion layer used
- For leisure and comfort shoes

FEATURES

- Lightweight natural cork mix
- Can be adapted to suit a wide range of foot proportions and shoe shapes
- Long-sole and short-sole blanks available
- Cork base also available with thermal reinforcement for extra stability
- High-quality sleek leather ensures better conditions within the shoe

INFORMATION

MODIFICATIONS

- Must be ground (foot orthosis shape at the sides and foot orthosis thickness)
- Easy to grind (grinding belt, grit 40–100, polishing brush on the blade disk)
- Easy to glue (standard workshop adhesives)

CLEANING / CARE

- Not washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

TOP COVER	CUSHION LAYER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 2 mm	–	✓	Wide, long, without reinforcement	20–32 ² 33–48	3 01650 5 00000**
Sleek leather, firm 2.25 mm	–	✓	Wide, long, thermal reinforcement	20–34 ² 35–48	3 01650 5 04721**
Sleek leather, firm 2.25 mm	–	Latex pad	Wide, long, thermal reinforcement	35–48	3 37650 5 22211**
Sleek leather, firm 2.25 mm	–	✓	Wide, long, wooden spline	35–48	3 37850 5 00000**
Sleek leather, soft 2 mm	TECmoll gray 1.5 mm (35° Shore A)	✓	Wide, long, thermal reinforcement	35–48 ¹	3 37050 5 82001**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing

GloboTec® Flexible Cork



ITEM NUMBER
3 03630 5 00000**

USES

- Cork / leather orthotic blank to provide stable support for the foot arches
- For the treatment of typical foot deformities, e.g. arch decreased / splay foot
- For leisure and comfort shoes

FEATURES

- Semi-flexible, lightweight foot orthosis base layer
- Can be adapted to suit a wide range of foot proportions and shoe shapes
- Long-sole and short-sole blanks available
- Cork base also available with thermal reinforcement for extra stability
- High-quality sleek leather ensures better conditions within the shoe

INFORMATION

MODIFICATIONS

- Must be ground (foot orthosis shape at the sides and foot orthosis thickness)
- Easy to grind (grinding belt, grit 40–100, polishing brush on the blade disk)
- Easy to glue (standard workshop adhesives)

CLEANING / CARE

- Not washable
- Do not expose to direct heat (heater, direct sunlight)
- We recommend that the functional properties are checked after six months

TOP COVER	METATARSAL PAD	VERSION	SIZE (**)	ITEM NUMBER
Sleek leather, soft 2 mm	Latex pad	Wide, long, thermal reinforcement	20–34 ¹⁺² 35–50	3 03630 5 22866**
Sleek leather, firm 2.25 mm	✓	Narrow, long, thermal reinforcement	35–46	3 03640 5 00000**
Sleek leather, firm 2.25 mm	✓	Medium width, long, thermal reinforcement	20–34 ² 35–48	3 03630 5 00000**
Sleek leather, firm 2.25 mm	✓	Wide, long, thermal reinforcement	20–34 ¹⁺² 35–50	3 03650 5 00000**
Sleek leather, firm 2.25 mm	✓	Extra wide, long, thermal reinforcement	35–48	3 03680 5 00000**
Sleek leather, firm 2.5 mm	✓ From size 36	Wide, long, thermal reinforcement, moderate heel cup	20–48	3 03750 5 00000**
Sleek leather, firm 2.5 mm	✓ From size 36	Wide, long, thermal reinforcement, high heel cup	20–48	3 03790 5 00000**

¹ Not a stock item. Please note that a longer delivery time will apply.

² Double sizing



VISCO

VISCOELASTIC FOOT ORTHOSES

FOR COMBATING ACUTE PROBLEMS

ViscoHeel®*

ViscoHeel® K*

Medical aid no. ViscoHeel:
08.03.06.0001

Supplied in pairs.



USES

- For the acute treatment of achillodynia and Haglund's deformity (ViscoHeel), varus and / or valgus misalignment of the heel (ViscoHeel K)
- Reduces shock loads on the ankle, knee, hip, and spine
- Soft cushioning reduces pressure symptoms under the heel
- Can also be used as pronation or supination cushion (ViscoHeel K)

FEATURES

- Viscoelastic, skin-friendly material
- Anatomically contoured
- Surface sloping at the sides (ViscoHeel K)

INFORMATION

CLEANING / CARE

- Can be cleaned with lukewarm water and a mild detergent
- Do not expose to direct heat (heater, direct sunlight)

ViscoHeel®

SIZE	SHOE SIZE	ITEM NUMBER
0	28–31	1 301702 2 000000
1	32–35	1 301702 2 000001
2	36–39	1 301702 2 000002
3	40–43	1 301702 2 000003
4	44–47	1 301702 2 000004


ViscoHeel® K

SIZE	SHOE SIZE	ITEM NUMBER
1	25–32	1 301710 2 000001
2	33–40	1 301710 2 000002
3	41–47	1 301710 2 000003

* Please refer to the relevant instructions for use for all information.

ViscoSpot^{®*}

Medical aid no.:
08.03.06.0008



USES

- For the acute treatment of heel spurs and achillodynia
- Targeted pressure relief in the area of the softest material (white spot)
- Additional relief of the insertion of the plantar fascia (blue spot)
- Recess along the course of the fascia

FEATURES

- Anatomical soft cushioning
- Smooth transition as a result of wavelike contours between the zones
- Reduces pressure at the edges
- Increased wearing comfort
- For use on either foot (with a compensating cushion for the other foot) or both feet

INFORMATION

CLEANING / CARE

- Can be cleaned with lukewarm water and a mild detergent
- Do not expose to direct heat (heater, direct sunlight)

ViscoSpot[®]

VERSION	SIZE	SHOE SIZE	ITEM NUMBER
Right / left	1	32–40	1 301720 2 000501
Right / left	2	41–47	1 301720 2 000502
Right	1	32–40	1 301720 2 000601
Right	2	41–47	1 301720 2 000602
Left	1	32–40	1 301720 2 000701
Left	2	41–47	1 301720 2 000702

* Please refer to the relevant instructions for use for all information.

ViscoBalance®*



USES

- For acute treatment of leg length differences and achillodynia
- Improves spine, knee, and hip joint symptoms caused by static strain
- Reduces shock loads on the ankle, knee, hip, and spine
- Suitable for leisure and sports shoes

FEATURES

- Viscoelastic, skin-friendly material
- Available in four sizes and three different heights: 3 mm, 5 mm and 10 mm
- Non-slip

INFORMATION

CLEANING / CARE

- Can be washed with lukewarm water and a mild detergent
- Do not expose to direct heat (heater, direct sunlight)

ViscoBalance®

SIZE	SHOE SIZE	ITEM NUMBER
1	32–35	1 30174* 0 050001
2	36–39	1 30174* 0 050002
3	40–43	1 30174* 0 050003
4	44–47	1 30174* 0 050004

*1 = height 3 mm, *2 = height 5 mm, *3 = height 10 mm

* Please refer to the relevant instructions for use for all information.

ViscoPed®*

ViscoPed® S*



ViscoPed®

Supplied in pairs.

USES

- Reduces load peaks
- Gentle correction of static or dynamic misalignment of the skeleton of the foot
- Special soft pads alleviate or eliminate painful irritation
- Reduces shock loads on the ankle, knee, hip, and spine
- Suitable for leisure and sports shoes

FEATURES

- Long-sole, viscoelastic foot orthosis
- Splay foot pad (ViscoPed) and slight longitudinal arch support
- Non-slip

INFORMATION

CLEANING / CARE

- Can be washed with lukewarm water and a mild detergent
- Do not expose to direct heat (heater, direct sunlight)

ViscoPed®

SIZE	SHOE SIZE	ITEM NUMBER
1	37–38	1 301730 1 000001
2	39–40	1 301730 1 000002
3	41–42	1 301730 1 000003
4	43–44	1 301730 1 000004
5	45–46	1 301730 1 000005

ViscoPed® S

SIZE	SHOE SIZE	ITEM NUMBER
0	35–36	1 381732 1 000000
1	37–38	1 381732 1 000001
2	39–40	1 381732 1 000002
3	41–42	1 381732 1 000003
4	43–44	1 381732 1 000004
5	45–46	1 381732 1 000005
6	47–48	1 381732 1 000006

* Please refer to the relevant instructions for use for all information.



MATERIAL

MATERIALS AND ADDITIONAL FEATURES

INDIVIDUALITY THANKS
TO CRAFTSMANSHIP

PedoPrint®

FOOTPRINT DEVICE*



ITEM NUMBER
3 60430 0 0000000

USES

- Image of the load-bearing area on the sole of the foot and the foot outline
- Allows for conclusions to be drawn about foot deformities
- Working basis for the provision of foot orthoses and shoes

FEATURES

- Traditional, portable footprint device for creating 2D footprints
- Available as a complete set including a rubber mat (with a grid pattern or plain), outline pen, ink roller, stamping ink (blue), 50 sheets of footprint paper, 10 sheets of protective film
- When the footprint device is used for a number of patients, the sheets of protective film ensure a better level of hygiene
- Individual components can also be ordered separately (see table)

INFORMATION

USE

- Use the ink roller to apply the ink thinly and evenly on one side of the rubber mat
- Put the sheet of footprint paper in place, attach the frame with the colored side of the rubber mat on the paper
- Position the foot on the clean side of the rubber mat and use the outline pen to draw round the foot at 90°
- Lift up the frame and remove the paper with the foot imprint
- Before taking a new imprint, ink must be reapplied to the rubber mat

CLEANING / CARE

- It is recommended that a surface disinfectant is used regularly to clean the footprint device and rubber mat

Complete PedoPrint® set

including: rubber mat (with a grid pattern or plain), outline pen, ink roller, stamping ink (blue), 50 sheets of footprint paper, 10 sheets of protective film

DESCRIPTION	VERSION	PACKAGING UNIT	ITEM NUMBER
Complete set	Rubber mat plain	Piece	3 60430 0 0000000
Complete set	Rubber mat with grid pattern	Piece	3 60430 0 0010000

PedoPrint® accessories

DESCRIPTION	VERSION	PACKAGING UNIT	ITEM NUMBER
Replacement rubber mat	Plain with plastic frame	Piece	3 60432 0 0000000
Replacement rubber mat	With grid pattern and plastic frame	Piece	3 60432 0 0010000
Ink roller	–	Piece	3 60431 0 0000000
Stamping ink	Blue	28 ml	3 60405 0 0000000
Outline pen	Synthetic material	Piece	3 60435 0 0000000
Sheets of footprint paper	Single-sided design	500 pieces	3 60440 0 0000000
Sheets of footprint paper	Double-sided design	500 pieces	3 60441 0 0000000
Sheets of protective film	–	100 pieces	3 60445 0 0000000

* Please refer to the relevant instructions for use for all information.

IMPRINT FOAM*



ITEM NUMBER
3 60450 1 0000000

USES

- Image of the load-bearing area on the sole of the foot and the foot outline
- Allows for conclusions to be drawn about foot deformities
- Working basis for the provision of foot orthoses and shoes

FEATURES

- Single-use PUR foam to create one pair of 3D foot imprints
- Available in a folding cardboard box, can be written on
- Imprint foam is available in the dimensions 325 × 145 × 38 mm and 325 × 145 × 50 mm (length × width × height)

INFORMATION

USE

- Open the imprint foam box
- If there are any wounds on the feet, cover them using sterile material
- With the patient either seated or standing, position one foot of the patient after the other on the foam and apply vertical pressure
- If necessary, mark the lateral line across the ball of the foot and the center of the ankle
- Remove the foot from the imprint foam and close the cardboard box
- To create foot orthosis lasts, plaster or last foam can be poured into the foam imprint




Imprint foam

DESCRIPTION	HEIGHT	PACKAGING UNIT	ITEM NUMBER
Imprint foam	38 mm	1 box with 50 units	3 60450 1 0000000
Imprint foam	38 mm	1 box with 25 units	3 60451 1 0000000
Imprint foam EXTRA	50 mm	1 box with 25 units	3 60457 1 0000000

 Manufacturer: Sprigron GmbH

* Please refer to the relevant instructions for use for all information.

METATARSAL PADS



FOAM RUBBER PADS

NO.	VERSION	SIZE (**)	ITEM NUMBER
1	Medium height, soft	1–4	3 60202 0 00*0000
2	High, soft	1–3	3 60207 0 00*0000
3	Flat, soft	1–2	3 60208 0 00*0000
4	Flat, soft	0–3	3 60209 0 00*0000
5	High, soft	1–4	3 60231 0 00*0000
6	High, very soft	1–5	3 60245 0 00*0000
7	Tear-shaped, leather cover, self-adhesive	1–4	3 60270 0 00*0000
8	T-shaped, leather cover, self-adhesive	1–4	3 60271 0 00*0000
9	High, soft	10–12	3 60650 2 0**0000
10	High, soft	1–3	3 60650 7 00*0000
11	Flat, soft, Anatomical shape	1–2	3 60650 3 00*0000
12	Medium height, soft, Anatomical shape	1–2	3 60650 4 00*0000

Packaging unit: 10 pairs each

CUSHIONS

FOREFOOT CUSHIONS | foam rubber | soft



LENGTH / Ø	SIZE	ITEM NUMBER
6 cm	1	3 60110 0 0010000
6.5 cm	2	3 60110 0 0020000
7 cm	3	3 60110 0 0030000
Packaging unit: 10 pairs each		

MIDFOOT CUSHIONS | foam rubber | retrocapital forefoot relief



LENGTH / Ø	SIZE	ITEM NUMBER
8 cm	10	3 60651 1 0100000
9 cm	11	3 60651 1 0110000
10 cm	12	3 60651 1 0120000
Packaging unit: 10 pairs each		

WEDGES

PES VALGUS WEDGES | with heel basket



LENGTH / Ø	SIZE	ITEM NUMBER
10 cm	0	3 60683 0 0000000
11 cm	1	3 60683 0 0010000
12 cm	2	3 60683 0 0020000
12.5 cm	3	3 60683 0 0030000

Packaging unit: 10 pairs each

ANGLED HEEL WEDGES | sponge rubber



LENGTH / Ø	SIZE	ITEM NUMBER
9.5 cm	1	3 60655 4 0010000°
7 cm	3	3 60655 4 0030000°

Packaging unit: 10 pairs each

SUPINATION WEDGE | EVA



LENGTH / Ø	SIZE	ITEM NUMBER
–	–	3 60310 0 0000900

Packaging unit: 10 pairs each

HEEL CUSHION

HEEL CUSHION | soft



LENGTH / Ø	SIZE	ITEM NUMBER
7.7 / 5.4 cm	1	3 60290 0 0010000
8.5 / 6.0 cm	2	3 60290 0 0020000
9.3 / 6.6 cm	3	3 60290 0 0030000
10 / 7.4 cm	4	3 60290 0 0040000

Packaging unit: 1 pair each

HEEL SPUR CUSHION

HEEL SPUR CUSHION | semi-firm | self-adhesive



LENGTH / Ø	SIZE	ITEM NUMBER
9.7 / 5.9 cm	10	3 60658 9 0100000
10.8 / 6.5 cm	12	3 60658 9 0120000

Packaging unit: 1 pair each

SMALL PADS

SMALL ADHESIVE PADS | beige | self-adhesive



LENGTH / Ø	ITEM NUMBER
33 mm	3 60016 0 0300500

Packaging unit: 50 pairs each

SMALL FLEECE PADS | beige | self-adhesive



LENGTH / Ø	ITEM NUMBER
33 mm	3 60017 0 0300500

Packaging unit: 50 pairs each

RUBBER FOOT ORTHOSES

SARDINIA



VERSION	SIZE (**)	ITEM NUMBER
Narrow	35–42	3 60280 2 **00000
Wide	38–46	3 60280 3 **00000

Packaging unit: 1 pair each

MALTA




VERSION	SIZE (**)	ITEM NUMBER
Narrow	36–42	3 60285 2 **00000
Wide	38–46	3 60285 3 **00000

Packaging unit: 1 pair each

TOP COVER MATERIALS

MICROFIBER



FEATURES

- Abrasion-resistant (Martindale tested in accordance with DIN EN ISO 12947)
- Comfortable, non-slip foot-side contact surface
- Polyester / polyethylene textile mix
- Selected products with silver ion coating
- Biocompatible in accordance with DIN EN ISO 10993-1

INFORMATION

MODIFICATIONS

- Easy to cut with shears or knife (non-directional)
- Easy to glue (standard workshop adhesives)

CLEANING / CARE

- Wipe using a mild cleanser and a sponge
- Air dry
- Do not expose to direct heat (heater, direct sunlight)


MICROFIBER

DESCRIPTION	THICKNESS	DIMENSIONS	ITEM NUMBER	
Microfiber blue with print (animal track pattern)	0.7 mm	1,000 × 1,400 mm	3 61223 5 0000 700	
Microfiber gray*	0.7 mm	1,000 × 1,400 mm	3 61225 0 0000 300	
Microfiber beige	0.7 mm	1,000 × 1,400 mm	3 61225 0 0000 500	
Microfiber brown*	0.7 mm	1,000 × 1,400 mm	3 61225 0 0000 600	
Microfiber blue*	0.7 mm	1,000 × 1,400 mm	3 61225 0 0000 700	
Microfiber navy*	0.7 mm	1,000 × 1,400 mm	3 61225 0 0001 800	
Microfiber anthracite*	0.7 mm	1,000 × 1,400 mm	3 61225 0 0002 200	
Microfiber macchiato*	0.7 mm	1,000 × 1,400 mm	3 61225 0 0002 600	

* with silver ion coating

TOP COVER MATERIALS

onSteam



FEATURES

- Abrasion resistant
- Comfortable foot-side contact surface
- onSteam available in a smooth and perforated design
- Biocompatible in accordance with DIN EN ISO 0993-1

INFORMATION

MODIFICATIONS

- Easy to cut with shears or knife (non-directional)
- Easy to glue (standard workshop adhesives)

CLEANING / CARE

- Wipe using a mild cleanser and a sponge
- Air dry
- Do not expose to direct heat (heater, direct sunlight)

onSteam

DESCRIPTION	THICKNESS	DIMENSIONS	ITEM NUMBER	
onSteam partly perforated, anthracite	0.8 mm	1,000 × 1,400 mm	3 61221 0 0002200	
onSteam partly perforated, black	0.8 mm	1,000 × 1,400 mm	3 61221 0 0000300	
onSteam partly perforated, beige	0.8 mm	1,000 × 1,400 mm	3 61221 0 0001500	
onSteam partly perforated, macchiato	0.8 mm	1,000 × 1,400 mm	3 61221 0 0000600	
onSteam partly perforated, navy	0.8 mm	1,000 × 1,400 mm	3 61221 0 0000700	

TOP COVER MATERIALS

LEATHER



FEATURES

- Abrasion resistant
- Comfortable foot-side contact surface
- Vegetable, chrome* tanned or combined** depending upon design
- Available in smooth and perforated design
- Available in various thicknesses

INFORMATION



MODIFICATIONS

- Easy to cut with shears or knife (non-directional)
- Easy to glue (standard workshop adhesives)
- Millable


CLEANING / CARE

- Air dry
- Do not expose to direct heat (heater, direct sunlight)



GOATSKIN

DESCRIPTION	THICKNESS	DIMENSIONS	ITEM NUMBER	
Goatskin	Approx. 0.8 mm	Whole hides	3 61674 0 0070 000	
Goatskin perforated, light, natural look, chrome tanned	Approx. 1 mm	Whole hides	3 61674 5 0070 000	

CALFSKIN LEATHER

DESCRIPTION	THICKNESS	DIMENSIONS	ITEM NUMBER	
Calfskin leather lining, natural look, perforated	0.8 – 1 mm	Whole hides	3 61619 0 0120 000	


SLEEK LEATHER HALVES

DESCRIPTION	THICKNESS	DIMENSIONS	ITEM NUMBER	
Sleek leather halves soft, perforated, natural look, vegetable tanned	Approx. 1.25 mm	Approx. 2.2 m ²	3 61611 0 0120 000	
Sleek leather halves firm, perforated, natural look, vegetable tanned	Approx. 1.25 mm	Approx. 2.2 m ²	3 61613 0 0120 000	

* No use of chromium VI salts! ** Chrome pretanned and then vegetable retanned

PLATE MATERIALS

TECfom / TECmoll / TEClast / TECstrong



USES

- Cushioning (TECfom, TECmoll) and construction (TECmoll, TEClast) materials for various types of foot orthosis
- For functional cushioning (TECfom, TECmoll), relief (TECfom) and support (TECmoll, TEClast)
- For sport, leisure and work

FEATURES

- Closed-cell soft foam (EVA / PE) (TECfom: 25° Shore A*; TECmoll: 35° Shore A*; TEClast: 45° Shore A*, TECstrong: 50° –55° Shore A*)
- Available in various thicknesses
- Biocompatible in accordance with DIN EN ISO 10993-1

INFORMATION

MODIFICATIONS

- Trimmable (sole leather shears / knife; saw can also be used for TECmoll and TEClast)
- Thermo-formable (TECfom / TECmoll: 80 to 120 °C hot air blower / oven; TEClast: 100 to 140 °C hot air blower / oven)
- Please note: Material shrinkage possible!
- Easy to grind (grinding belt: grit 24–100, pumice stone, blade disk)
- Easy to glue (first roughen extensively, standard workshop adhesives)

STORAGE





- Store hanging or lying
- Store at 10–25 °C
- Store in a dry place, do not expose to direct heat (heater, direct sunlight)

TECfom (EVA / PE)

DESCRIPTION	SHORE HARDNESS	THICKNESS	DIMENSIONS	ITEM NUMBER	
TECfom sapphire blue	25° Shore A	3 mm	1,150 × 1,150 mm	3 61457 0 0300 700	
TECfom light azure blue	25° Shore A	3 mm	1,150 × 1,150 mm	3 61467 1 0300 700	
TECfom sapphire blue	25° Shore A	6 mm	1,150 × 1,150 mm	3 61457 0 0600 700	
TECfom light azure blue	25° Shore A	6 mm	1,150 × 1,150 mm	3 61467 1 0600 700	
TECfom gray	25° Shore A	6 mm	1,150 × 1,150 mm	3 61467 1 0600 400	

* Production-related Shore hardness tolerance: +/-10%

TECmoll (EVA / PE)

DESCRIPTION	SHORE HARDNESS	THICKNESS	DIMENSIONS	ITEM NUMBER	
TECmoll silver gray	35° Shore A	3 mm	1,000 × 1,000 mm	3 61458 0 0300 400	
TECmoll silver gray	35° Shore A	6 mm	1,000 × 1,000 mm	3 61458 0 0600 400	
TECmoll silver gray	35° Shore A	15 mm	1,000 × 1,000 mm	3 61458 0 1500 400	
TECmoll sapphire blue	35° Shore A	6 mm	1,000 × 1,000 mm	3 61458 0 0600 700	

TEClast (EVA / PE)

DESCRIPTION	SHORE HARDNESS	THICKNESS	DIMENSIONS	ITEM NUMBER	
TEClast gray	45° Shore A	6 mm	1,000 × 1,000 mm	3 61459 0 0600 400	
TEClast gray	45° Shore A	15 mm	1,000 × 1,000 mm	3 61459 0 1500 400	

TECstrong (EVA / PE)

DESCRIPTION	SHORE HARDNESS	THICKNESS	DIMENSIONS	ITEM NUMBER	
TECstrong brown	55° Shore A*	6 mm	1,260 × 980 mm	3 614540 0600 600	

PLATE MATERIALS

TECsoft (PE)



USES

- Cushioning and construction material for various types of foot orthosis
- For functional cushioning, relief and support
- For leisure and work

FEATURES

- Closed-cell soft foam (PE, 25° Shore A*)
- Available in various thicknesses

INFORMATION

MODIFICATIONS

- Trimmable (sole leather shears / shears)
- Thermoplastic molding possible (100 °C hot air blower / oven)
- Please note: Material shrinkage possible!
- Easy to grind (grinding belt: grit 24–100, pumice stone, blade disk)
- Easy to glue (first roughen extensively, standard workshop adhesives)

STORAGE

- Store hanging or lying
- Store at 10–25 °C
- Store in a dry place, do not expose to direct heat (heater, direct sunlight)

TECsoft (PE)

DESCRIPTION	SHORE HARDNESS	THICKNESS	DIMENSIONS	ITEM NUMBER	
TECsoft gray	25° Shore A*	3 mm	1,200 × 1,000 mm	3 61440 1 0300 300	
TECsoft gray	25° Shore A*	5 mm	1,200 × 1,000 mm	3 61440 1 0500 300	

* Production-related Shore hardness tolerance: +/-10%

PLATE MATERIALS

TECtop (PE)



USES

- Cushioning and construction material for various types of foot orthosis
- For functional cushioning and support
- For sport, leisure and work

FEATURES

- Closed-cell soft foam (EVA / PE, 30° Shore A*)
- Available in various thicknesses

INFORMATION

MODIFICATIONS

- Trimmable (sole leather shears / shears)
- Thermoplastic molding possible (100–120 °C hot air blower / oven)
- Please note: Material shrinkage possible!
- Easy to grind (grinding belt: grit 24–100, pumice stone, blade disk)
- Easy to glue (first roughen extensively, standard workshop adhesives)

STORAGE

- Store hanging or lying
- Store at 10–25 °C
- Store in a dry place, do not expose to direct heat (heater, direct sunlight)

TECtop (PE)

DESCRIPTION	SHORE HARDNESS	THICKNESS	DIMENSIONS	ITEM NUMBER	
TECtop silver-gray, perforated	30° Shore A*	3 mm	1,100 × 1,100 mm	3 61453 0 0300 400	
TECtop red, perforated	30° Shore A*	4 mm	1,100 × 1,100 mm	3 61453 0 0400 800	
TECtop earth brown, perforated	30° Shore A*	4 mm	1,100 × 1,100 mm	3 61453 0 0400 600	
TECtop sapphire blue, perforated	30° Shore A*	2 mm	1,100 × 1,100 mm	3 61453 0 0200 700	
TECtop sapphire blue, perforated	30° Shore A*	4 mm	1,100 × 1,100 mm	3 61453 0 0400 700	

* Production-related Shore hardness tolerance: + / -10%

PLATE MATERIALS

PORON / GLOBOFIT (PU)



USES

- Cushioning material for various types of foot orthosis
- For functional cushioning and relief
- For sport, leisure and work

FEATURES

- Open-cell soft foam (PU)
- Breathable
- Microcellular structure
- Fungus-resistant
- Available in various thicknesses
- Biocompatible in accordance with DIN EN ISO 10993-1

INFORMATION

MODIFICATIONS

- Easy to cut with shears
- Easy to grind (grinding belt: grit 24–100, pumice stone, blade disk)
- Easy to glue (first roughen extensively, standard workshop adhesives)




STORAGE

- Store hanging or lying
- Store at 10–25 °C
- Store in a dry place, do not expose to direct heat (heater, direct sunlight)

PORON

DESCRIPTION	THICKNESS	DIMENSIONS	ITEM NUMBER	
Poron 2001 salmon pink, roughened on one side, very soft	3 mm	1,370 × 1,000 mm	3 61445 0 0300 100	

GLOBOFIT

DESCRIPTION	SHORE HARDNESS	THICKNESS	DIMENSIONS	ITEM NUMBER	
Globofit beige, roughened on both sides	Approx. 12° Shore A*	3 mm	1,370 × 1,000 mm	3 61443 0 2300 500	
Globofit light blue, roughened on both sides	Approx. 12° Shore A*	3 mm	1,370 × 1,000 mm	3 61443 0 2300 700	
Globofit light blue, roughened on both sides	Approx. 12° Shore A*	6 mm	1,370 × 1,000 mm	3 61443 0 2600 700	

* Production-related Shore hardness tolerance: +/-10%

PLATE MATERIALS

TECesd (EVA / PE)



USES

- Construction material for various types of foot orthosis
- For functional support
- For sport, leisure and work

FEATURES

- Closed-cell soft foam (EVA / PE, 40° Shore A*)
- Available in various thicknesses
- High conductivity

INFORMATION

MODIFICATIONS

- Trimmable (sole leather shears / shears / saw)
- Thermoplastic molding possible (100–120 °C hot air blower / oven)
- Please note: Material shrinkage possible!
- Easy to grind (grinding belt: grit 24–100, pumice stone, blade disk)
- Easy to glue (first roughen extensively, standard workshop adhesives)

STORAGE

- Store hanging or lying
- Store at 10–25 °C
- Store in a dry place, do not expose to direct heat (heater, direct sunlight)

TECesd (EVA / PE)

DESCRIPTION	SHORE HARDNESS	THICKNESS	DIMENSIONS	ITEM NUMBER	
TECesd black	40° Shore A	6 mm	1,050 × 1,050 mm	3 61106 0 0600 300	
TECesd black	40° Shore A	10 mm	1,050 × 1,050 mm	3 61106 0 1000 300	

* Production-related Shore hardness tolerance: + / -10%

ADDITIONAL BAUERFEIND PRODUCT CATEGORIES



VenoTrain COMPRESSION STOCKINGS

are characterized by clinically proven effectiveness and unique material properties. They improve venous activity and counteract venous leg problems. For both mild to moderate and advanced vein problems, as well as for work or travel, the VenoTrain product range offers a suitable stocking. They can easily be combined with various styles, thanks to a wide range of colors and designs.



COMPRESSION STOCKINGS FOR EDEMA TREATMENT

Every day is more pleasant when you feel good. That is why the flat knit, custom-made VenoTrain® curaflow compression stockings were developed. This versatile and colorful product perfectly adapts to the individual shape of the lymphedema or lipedema patient's legs, arms, hands, or thorax. VenoTrain® curaflow provides noticeable relief so patients can live their everyday lives feeling carefree, confident, and attractive.



BAUERFEIND'S SUPPORTS AND ORTHOSES

help after injuries as well as with acute and chronic orthopedic conditions.

They improve the feeling of stability and provide confidence during everyday activities and exercise. The supports and orthoses are made from lightweight and breathable materials to perfectly underpin therapeutic exercise.

They are easy to handle during everyday use, thanks to innovative solutions, such as quick-fit fasteners and integrated donning aids, and will win you over with their premium quality.

MEDICAL AIDS
HAVE BEEN
PROVEN TO BE
EFFECTIVE*

*** Supports / orthoses, compression stockings and foot orthoses help patients get back on their feet quickly after an injury or in case of disorders.**

(Source: representative patient survey conducted by the Allensbach Institute for Public Opinion Research on behalf of eurocom e.V., 2023, n=n=1,300 users of supports / orthoses, compression stockings, and foot orthoses)

For an overview of all treatment areas, please visit:
experts.bauerfeind.com

BAUERFEIND AG TERMS AND CONDITIONS OF SALE AND DELIVERY (GTC)

VALID FOR: GERMANY

Prices

The prices quoted are net prices in Euros ex works and do not include statutory VAT.

Shipping

Shipping is at the risk and on the account of the purchaser, even if the purchaser does not bear the costs of transportation. Business interruptions affecting us or any of our suppliers, lack of suitable personnel, strikes, lockouts and force majeure are causes for extending delivery times and relieve us of any delivery obligation for the duration of the hindrance.

Shipping costs

We calculate our shipping costs as a flat rate of EUR 7.00 for a net total order value of up to EUR 50.00, and EUR 4.00 for a net total order value of up to EUR 200.00. No shipping costs need to be paid for custom-made solutions. Special shipping requirements will be charged at cost.

Returns

Returns are only accepted outside the warranty terms if we have been notified in advance and approved the return (notification) with inclusion of the shipping documents and under the assumption that the returned goods are in faultless condition. We reserve the right to reject any returns.

Except where covered by warranty, we charge for processing returns. For returns within one month of the delivery date, this charge is EUR 10.00 for each return. For goods returned one month or more after the delivery date, a handling fee of 20 percent of the net order value, or a maximum of EUR 25.00 per product and unit, is charged; for returns six to twelve months after the delivery date, the fee is 40 percent of the net order value, or a maximum fee of EUR 50.00 for each product and unit.

Except where covered by warranty, no returns, including, but not limited to, custom-made orders, will be accepted after twelve months following the delivery date.

Payments

Payments are due immediately following delivery. We offer a 4 percent discount for payments made by direct debit and a 3 percent discount for payments made within ten days of the invoice date. Otherwise, the payment terms are 30 days net. The buyer shall be considered in arrears once 30 days have passed from the invoice date. From this point on, any further deliveries will only be made against advance payment.

From the date the customer is in arrears, we will charge default interest of 8 percent above the applicable base interest rate. Payments by check, bank draft, etc. are deemed to be made only when the amount has been credited fully and irrevocably to our account.

Retention of ownership

We retain the proprietary rights to all of the delivered goods until the purchaser has settled all of the outstanding accounts due to us as a result of the business relationship. The retention of title shall also apply pro rata if the goods in question have been processed or modified. If the goods are resold, all third-party claims arising therefrom shall be transferred to us without reservation up to the total amount owed to us. Notwithstanding the assignment, the buyer is entitled to the collection of these amounts. Our right to collect the amounts is not affected by the preceding provision; however, we undertake to refrain from collection unless the buyer defaults on payment. If the buyer defaults on payment, we can demand that the buyer disclose the assigned claims and the names of the debtors, provide all information required for collection of the debt, deliver all required documents and inform the debtors (third party) of the assignment.

The buyer shall not be entitled to pledge or assign the goods as security. The buyer shall notify us without undue delay about any seizures, confiscations, or any other dispositions taken by a third party, and shall make available all information and documents necessary for us to defend our rights.

We undertake to release the security owed to us at the buyer's request to the extent that its value does not exceed that of the claim to be secured by more than 20 percent and provided the claim in question has not yet been settled.

Complaints

Warranty claims in respect of obvious defects shall be void if we do not receive immediate written notification of these claims in accordance with Section 377 of the German Commercial Code (HGB). In all other cases, the period for warranty claims is one year after handing over the purchased goods. If the goods are defective, we may replace or repair the goods at our discretion. For reasons of hygiene, repairs will only be performed on products that have been washed or are completely hygienic. Upon receipt of unwashed goods, we reserve the right to claim a suitable sum for additional work. Provided that we fulfill our

supplementary performance obligations, especially in regard to rectifying defects or delivering a defect-free product, the customer has no right to demand a price reduction or to withdraw from the contract, unless we fail to deliver our supplementary performance. Defects which are caused by normal wear and tear, overuse, misuse, or neglect of care are not covered by warranty.

Liability

We are liable only in cases of gross negligence and intent as well as for breaching obligations that make the proper performance of the contract possible in the first place and on the compliance with which the buyer may regularly depend ("material obligations").

Regarding the negligent violation of a material obligation, our liability is limited to damages that are foreseeable when the contract is signed and that are typical of this kind of contract.

We are not liable for the negligent violation of contractual obligations which are not material obligations.

Provided that our liability is limited or excluded, this shall also apply to the liability of our employees, representatives, or vicarious agents.

The liability limitations or exclusions above do not apply in cases of fraudulent concealment of defects, assumption of a guarantee or procurement risk, for liability based on the Product Liability Law, as well as for physical injuries (injury to life, limb, or health). No change of burden of proof to the disadvantage of the buyer is associated with this provision.

Applicable law and place of jurisdiction

These Terms and Conditions shall be governed by the laws of the Federal Republic of Germany, excluding the United Nations Convention on Contracts for the International Sale of Goods (CISG).

The place of jurisdiction for all disputes resulting from or associated with this contract shall be Gera, provided that the buyer is a legal entity under public law or a special fund under public law within the meaning of the German Commercial Code. However, we are also entitled in all cases to file suit at the place of fulfillment or at the general jurisdiction of the user. Overriding legal regulations, particularly regarding exclusive jurisdictions, remain unaffected.

The place of performance for the delivery is our company headquarters.

Severability clause

If a provision in these Terms and Conditions is or becomes fully or partly void, ineffective or unenforceable ("Flawed Provision"), the effectiveness and the enforceability of the remaining provisions in these Terms and Conditions will not be affected. The parties agree at this time to find a provision instead of the flawed provision that, within the limits of what is legally possible, is as close to what the parties would have agreed for the meaning and purpose of this provision if they had known that the provision was flawed. If the inaccuracy of the provision is based on a measure of performance or time (period or deadline) specified therein, the provision must be agreed with a legally permissible measure that most closely reflects the original measure. The same applies to potential loopholes in these Terms and Conditions. This severability clause does not entail the mere reverse onus, but also waives Section 139 of the German Civil Code in general.

Written form

Changes and/or additions to these Terms and Conditions must be made in writing to take effect (handwritten signature); this shall also apply to any changes and additions to this written form requirement.

Data protection clause / conversation recording

For quality purposes, we randomly record calls between the buyer and our employees. The phone recordings solely serve the purpose of improving our service and are deleted after evaluation. The data will not be disclosed to third parties. The buyer hereby consents to the recording of telephone calls for that sole purpose. Consent is voluntary and may be revoked by the buyer at any time with effect for the future.

Please see www.bauerfeind.de/ds for further information on the latest version of Bauerfeind AG's General Data Protection Policy.

A BETTER QUALITY OF LIFE FOR YOUR CUSTOMERS!

Disclaimer

All information in this brochure is up to date at the time of printing. This brochure – depending on its version – gives an initial insight into this area of expertise, and suggests selected medical diagnosis and treatment options. Even though the details are scientifically substantiated, we do not claim that all information is complete or that all schools of thought are represented. All contents in this brochure or its online version, including extracts, in particular those about diagnostics and treatment options (texts, images, graphs, etc.), are understood as non-binding support for the treatment that must be chosen by a physician or professional. This brochure is no substitute for the necessary personal medical history and its overall assessment by a physician/professional, considering ongoing treatment, if applicable. The contents in this brochure are exclusively designed for professionals, not for anyone unrelated to the subject, in particular not for patients. If in doubt, please seek additional advice from a specialist. We are also happy to help.

© Bauerfeind AG/2023

Do you have any questions, suggestions or requests?

Please contact us:

BAUERFEIND AG

Triebeser Straße 16
07937 Zeulenroda-Triebes
Germany

T +49 (0) 36628 66-35 00

F +49 (0) 36628 66-39 99

E info@bauerfeind.com

→ [BAUERFEIND.COM](https://www.bauerfeind.com)